



Diagnostics
is in our blood.



Catalogue

Educational tools

Digital

Webinars

An innovative interactive learning tool reinforcing Stago's commitment to medical education and training in the field of Thrombosis and hemostasis.

- 4 times a year since 2014
- With the contribution of world-renowned experts


Sign up to be informed & participate in each «live» webinar and to access previous sessions «on-demand».

Podcasts

Ask Stago is the Stago podcast series at destination of the Laboratory staff around the daily work at the hemostasis bench.

Every month and in less than 10 min, our team wants to provide "expert answers to your expert questions". Content is scientific and technical in nature. It is intended as an educational tool for laboratory professionals and topics discussed are not intended as recommendations or as commentary on appropriate clinical practice. Available on Stago youtube channel and all podcast platforms

www.stagowebinars.com



iHemostasis (for tablet only)



A global educational tool dedicated to the pathophysiology of blood coagulation illustrated by clinical cases, designed for clinical pathologists, physicians, medicine and pharmacy students, clinical laboratory staff, patient care personnel or anyone wishing to learn more in the field of hemostasis and thrombosis.

5 INTERACTIVE SECTIONS TO LEARN

- Coagulation cascade
- Clinical cases
- Practical Manual series
- Special Focus
- Quick guide



Haemoscore (for smartphone & tablet)



An application developed with an international expert panel.

- Compiles, in a clear and simple way, the most recognized and useful **clinical scores and diagnosis algorithms** in the field of thrombosis and hemostasis.
- **42 scores and algorithms related to:** DVT & PE, pregnancy, atrial fibrillation, bleeding, antithrombotic therapy and other thrombotic conditions.

iHemostasis



For tablets only

Haemoscore



Tablet & smartphone

Ask Stago



Available on podcast platforms

Digital links

- Available 24/7: you can download reagent package inserts and instrumentation user documents with Intended use, Test principle, Clinical application and Performance.
- Available in more than 20 languages*
- Document traceability with historical versions of each package insert at hand when you need it.
- Assistance in meeting accreditation requirements

Instruction For Use



Available on Stago.com

Webinars



Available on Stagowebinars.com

Haemostasis



General information on Stago.com

*Availability depending upon the country.

SUMMARY

SYSTEMS

STA R Max systems	4
STA Workcell Max	8
STA Compact Max systems	9
STA Satellite Max	12
STart Max	14
ST Genesisia system	16
Calibrated Automated Thrombogram	18
Platelet Aggregometers TA-8V & TA-4V	19

DIGITAL SOLUTIONS

Coag.One by Stago	22
Connect.One by Stago	24
My Expert QC	26
My Qualiris QC	28

REAGENTS

SCREENING TESTS	31	THROMBOPHILIA	46
Prothrombin Time (PT)	32	Major Physiological Inhibitors	46
Activated Partial Thromboplastin Time (APTT)	32	Antiphospholipid syndrome	47
Fibrinogen	33	Other Inhibitors	47
Thrombin Time (TT)	33	PRIMARY HEMOSTASIS	49
FACTOR ASSAYS	36	Von Willebrand Disease	49
Deficient plasmas	36	Platelets	50
Other methods	38	Monoclonal Antibodies	52
FIBRIN FORMATION & DEGRADATION	40	GLOBAL TESTS	54
Fibrin formation	40	Thrombin Generation	54
Degradation products	40	CALIBRATORS & CONTROLS	55
FIBRINOLYSIS	42	Calibration plasmas	55
Fibrinolysis	42	Control plasmas	56
ANTICOAGULANT TREATMENTS	43	Correlation Plasma sets	57
Heparins & Direct Anti-Xa	43	External Quality Control Program	58
HIT	44	OTHER REAGENTS	59
DTI	44	Auxiliary Reagents	59

RESEARCH PRODUCTS

Research kits	62		
Purified proteins	64		
Chromogenic substrates & Activators	64		
Monoclonal Antibodies for Flow Cytometry	65		
Thrombin Generation (CAT)	66		
Thrombin Generation (ST Genesisia)	66		
ExpertCor	66		
Heparin Induced Thrombocytopenia	66		
		Index	68

SYSTEMS



This catalogue contains information on products which are targeted to a wide range of audiences and could contain product details or information concerning goods which are not accessible or valid in your country.



STA R MAX SYSTEMS

Third version of the high-end instrument from the "Max Generation" line, the STA R Max 3 is a fully automated analyzer which delivers sample integrity verification to ensure quality of results, thanks to the EPC module and Viscosity-Based Detection System.

Supplied with

- 1 accessory pack
- 1 21.5" LCD touch screen
- 1 computer table
- STA Coag Expert (optional)
- Visual Alarm Device (on STA R Max 3 only - optional on STA R Max)

To order*

120 V 60 Hz

- 58413** STA R Max
- 58412** STA R Max cap piercing
- 58411** STA R Max & STA Coag Expert
- 58410** STA R Max cap piercing & STA Coag Expert

230 V 50 Hz

- 58407** STA R Max
- 58406** STA R Max cap piercing
- 58405** STA R Max & STA Coag Expert
- 58414** STA R Max cap piercing & STA Coag Expert

- 58404** STA R Max cap piercing & STA Coag Expert
- 59304** STA R Max 3
- 59305** STA R Max 3 cap piercing
- 59306** STA R Max 3 & STA Coag Expert
- 59307** STA R Max 3 cap piercing & STA Coag Expert

230 V 60 Hz

- 58417** STA R Max
- 58416** STA R Max cap piercing
- 58415** STA R Max & STA Coag Expert

*Availability depending upon the country.

Parameters

- PT (Quick and Owren PT)
- APTT
- Fibrinogen
- Thrombin Time
- Reptilase Time
- Extrinsic pathway factors
- Intrinsic pathway factors
- Factor XIII
- Anti-Xa (UFH, LMWH, Rivaroxaban...)
- Anti-IIa
- D-Dimers, Fibrin Monomers and Fibrin Degradation Products
- Antithrombin
- Protein C
- Activated Protein C Resistance
- Protein S
- Lupus Anticoagulants
- VWF Ag & RCo
- Microparticles
- Plasminogen, Antiplasmin and TAFI
- Calibrators
- Quality Controls

Olivier Q. - Middleware Team Manager

"The health ecosystem is rapidly evolving and the management of the data generated by the IVD device and software is more and more crucial. To proactively support our customers, Stago has integrated new skills and positions such as the middleware specialists team created in 2012."



STA R MAX LINE / STA R / STA R EVOLUTION EXPERT SERIES Solutions & Disposables

Cat. NR.	Product name	Product description	Pack.
Solutions			
00975	STA-Desorb U	24 vials of STA-Desorb U	24 x 15 mL
00973	STA-Cleaner Solution	6 vials of STA-Cleaner Solution	6 x 2.5 L
Disposables			
38669	STA-Cuvettes	6 rolls of 1000 STA-Cuvettes	6 x 1000
00802	STA-Microcups Microtubes in siliconed glass for reagents control and calibrator plasmas	1 box of 100 STA-Microcups	1 x 100
00741	STA - Microtainer Plastic cups for samples	1 box of 500 Microcontainer	1 x 500
00797	STA-Mini Reducer Plastic chimneys for reagents (from 4 to 6 mL vials) - Improved stability	1 box of 100 STA-Mini Reducer	1 x 100
00801	STA-Maxi Reducer Plastic chimneys for reagents (from 8 to 15 mL vials) - Improved stability	1 box of 100 STA-Maxi Reducer	1 x 100
00760	STA-Micro Reducer Plastic chimneys for reagents (from 1 to 2 mL vials) - Improved stability	1 box of 100 STA-Micro Reducer	1 x 100
27425	White stirring magnet (for STA-Néoplastine and STA-NeoPTimal)	1 item	1
26674	Red stirring magnet (for STA-C.K. Prest)	1 item	1
89174	White stirring magnet (for STA-Néoplastine and STA-NeoPTimal)	10 items	1 x 10
89238	Red stirring magnet (for STA-C.K. Prest)	10 items	1 x 10

*Availability depending upon the country.

STA R MAX LINE / STA R / STA R EVOLUTION EXPERT SERIES

Disposables spare parts

Cat. NR.	Product name	Product description	Pack.
89156	Needle 1 Quick change for cap piercing on STA R Max 2 & 3	1 item	1
39164	Needle 1 for cap piercing	1 item	1
39249	Needle 1	1 item	1
39250	Needle 2	1 item	1
27307	Needle 3	1 item	1
27530	Teflon tips for syringe and o-ring	1 item	2 x 6
26699	Halogen lamp	1 item	1
27458	Liquid filter	1 item	1
39738	Liquid filter	10 items	1 x 10
27538	Hamilton syringe and o-ring	1 item	1
87063	Suction head	2 items	1 x 2
38517	Air filter	1 item	1
39880	5 x 20 T 1A fuses	10 fuses	1 x 10
26682	5 x 20 T 2A fuses	10 fuses	1 x 10
26684	5 x 20 T 5A fuses	10 fuses	1 x 10

STA R MAX LINE / STA R / STA R EVOLUTION EXPERT SERIES

Disposables spare parts

Cat. NR.	Product name	Product description	Pack.
26681	5 x 20 T 6.3A fuses	10 fuses	1 x 10
39863	5 x 20 T 10A fuses	10 fuses	1 x 10
27576	6.3 x 32 T 2A fuses	10 fuses	1 x 10
27578	6.3 x 32 T 5A fuses	10 fuses	1 x 10
27575	6.3 x 32 T 6.25A fuses	10 fuses	1 x 10
27224	6.3 x 32 T 10A fuses	10 fuses	1 x 10
26694	6.3 x 32 T 15A fuses	10 fuses	1 x 10
38125	Colorimetric box filter	1 item	1
39134	Wheel for tube rotation	2 items	1 x 2
80343	Soft stylet kit	10 items	1 x 10

*Availability depending upon the country.

STA R MAX LINE / STA R / STA R EVOLUTION EXPERT SERIES

Auxiliary materials

Cat. NR.	Product name	Product description	Pack.
26555	Ball extractor	1 item	1
27543	Microcups reducer	2 items	1 x 2
39002*	Aluminium microcontainer reducer (microvolume)	10 items	1 x 10
39010*	Sample tray	2 items	1 x 2
89086**	Sample tray	2 items	1 x 2
38777	Sample rack holder (for 5 racks)	1 item	1
39772	Sample rack holder (for 5 racks)	10 items	1 x 10
39147	Sample rack	10 items	1 x 10
89381	Sample rack 11-13	10 items	1 x 10
89257	Microcontainer standard rack	2 items	1 x 2
89318	Screw microtube kit	2 screw microtube racks • 10 microtube adaptors • 1 mapping tool	1 x 2 1 x 10
89377	Visual Alarm Device	1 item	1



*For STA-R Evolution only

**For STA-R Max only

*Availability depending upon the country.



STA WORKCELL MAX

The STA Workcell Max is a transformative specimen processing platform designed to meet the rigorous requirements of efficient hemostasis testing. By automating front-end specimen processing, rack preparation and specimen routing, the STA Workcell Max enables allocation of skilled laboratory technologists for other more complex, value-added tasks.

To order*

- 87506 Loading Module
- 87510 Hettich Centrifuge
- 87507 Dispatcher Module
- 87513 Dispatcher Module Stand Alone
- 87508 Long Track Module
- 87511 Short Track Module
- 87509 Instrument Transfer Module
- 87512 Output Module



*Availability depending upon the country.



STA COMPACT MAX SYSTEMS

Dedicated for mid-volume throughput laboratories, the STA Compact Max 3 is a fully automated analyzer which delivers sample integrity verification to ensure quality of results, thanks to the EPC module and Viscosity-Based Detection System.

Supplied with

- 1 accessory pack
- 1 21.5" LCD touch screen
- STA Coag Expert (optional)
- Visual Alarm Device (optional)

To order*

- 58712** STA Compact Max
- 58713** STA Compact Max cap piercing
- 58732** STA Compact Max with STA Coag Expert
- 58733** STA Compact Max cap piercing with STA Coag Expert
- 59202** STA Compact Max 3
- 59203** STA Compact Max 3 cap piercing
- 59206** STA Compact Max 3 with STA Coag Expert
- 59207** STA Compact Max 3 cap piercing with STA Coag Expert

Parameters

- PT (Quick and Owren PT*)
- APTT
- Fibrinogen
- Thrombin Time
- Reptilase Time
- Extrinsic pathway factors
- Intrinsic pathway factors
- Factor XIII
- Anti-Xa (UFH, LMWH, Rivaroxaban...)
- Anti-IIa
- D-Dimers, Fibrin Monomers and Fibrin Degradation Products
- Antithrombin
- Protein C
- Activated Protein C Resistance
- Protein S
- Lupus Anticoagulants
- VWF Ag
- Microparticles
- Plasminogen, Antiplasmin and TAFI
- Calibrators
- Quality controls

*Availability depending upon the country.

STA COMPACT MAX LINE / STA COMPACT Solutions & Disposables

Cat. NR.	Product name	Product description	Pack.
Solutions			
00975	STA-Desorb U	24 vials of STA-Desorb U	24 x 15 mL
00973	STA-Cleaner Solution	6 vials of STA-Cleaner Solution	6 x 2.5 L
Disposables			
38669	STA-Cuvettes	6 rolls of 1000 STA-Cuvettes	6 x 1000
00802	STA-Microcups Microtubes in siliconed glass for reagents control and calibrator plasmas	1 box of 100 STA-Microcups	1 x 100
00741	STA - Microtainer Plastic cups for samples	1 box of 500 Microcontainer	1 x 500
00797	STA-Mini Reducer Plastic chimneys for reagents (from 4 to 6 mL vials) - Improved stability	1 box of 100 STA-Mini Reducer	1 x 100
00760	STA-Micro Reducer Plastic chimneys for reagents (from 1 to 2 mL vials) - Improved stability	1 box of 100 STA-Micro Reducer	1 x 100
00801	STA-Maxi Reducer Plastic chimneys for reagents (from 8 to 15 mL vials) - Improved stability	1 box of 100 STA-Maxi Reducer	1 x 100
27425	White stirring magnet (for STA-Néoplastine and STA-NeoPTimal)	1 item	1
26674	Red stirring magnet (for STA-C.K. Prest)	1 item	1
89174	White stirring magnet (for STA-Néoplastine and STA-NeoPTimal)	10 items	1 x 10
89238	Red stirring magnet (for STA-C.K. Prest)	10 items	1 x 10

*Availability depending upon the country.

STA COMPACT MAX LINE / STA COMPACT

Disposables spare parts

Cat. NR.	Product name	Product description	Pack.
89156	Needle 1 Quick change for cap piercing on STA Compact Max 3	1 item	1
39022	Needle 1 for cap piercing	1 item	1
38646	Needle 1	1 item	1
27354	Needle 2	1 item	1
27307	Needle 3	1 item	1
27530	Teflon tips for syringe and o-ring	1 item	2 x 6
27458	Liquid filter	1 item	1
39738	Liquid filter	10 items	1 x 10
26699	Halogen lamp	1 item	1
27538	Hamilton syringe and o-ring	1 item	1
27421	Suction head	2 items	1 x 2
27420	Air filter	2 items	1 x 2
80675	Rising plate assay air filter	1 item	1
26538	Colorimetric box filter	1 item	1
38125	Colorimetric box filter	1 item	1
87018	Tubing EV piercing	6 items	1 x 6
80343	Soft stylet kit	10 items	1 x 10
27034	6.3 x 32 T 0.3A fuses	10 fuses	1 x 10
27579	6.3 x 32 T 1A fuses	10 fuses	1 x 10
27037	6.3 x 32 T 4A fuses	10 fuses	1 x 10
27575	6.3 x 32 T 6.25A fuses	10 fuses	1 x 10
27225	6.3 x 32 T 8A fuses	10 fuses	1 x 10
27224	6.3 x 32 T 10A fuses	10 fuses	1 x 10
39783	5 x 20 T 315mA fuses	10 items	1 x 10
39880	5 x 20 T 1A fuses	10 fuses	1 x 10
39784	5 x 20 T 4A fuses	10 items	1 x 10
39785	5 x 20 T 8A fuses	10 items	1 x 10
39863	5 x 20 T 10A fuses	10 fuses	1 x 10
26681	5 x 20 T 6.3A fuses	10 fuses	1 x 10

STA COMPACT MAX LINE / STA COMPACT

Auxiliary materials

Cat. NR.	Product name	Product description	Pack.
26555	Ball extractor	1 item	1
39002	Aluminium microcontainer reducer (microvolume)	10 items	1 x 10
27543	Microcups reducer	2 items	1 x 2
88618	QWERTY keyboard for STA Compact Max	1 item	1
89377	Visual Alarm Device	1 item	1

*Availability depending upon the country.



STA SATELLITE MAX

The STA Satellite Max is dedicated for low-volume throughput laboratories. Its design has been developed for continuous result reporting and minimal user maintenance - using Viscosity-Based Detection System.

Supplied with

- 1 accessory pack
- Screen Touch LCD colour screen 17" + External PC

To order*

58030 STA Satellite Max

Parameters*

- PT (Quick and Owren PT)
- APTT
- Fibrinogen
- Thrombin Time
- Anti-Xa (UFH, LMWH)
- D-Dimers
- Antithrombin
- Calibrators
- Quality controls

*Availability depending upon the country.

STA SATELLITE MAX / STA SATELLITE USB Solutions & Disposables

Cat. NR.	Product name	Product description	Pack.
Solutions			
00975	STA-Desorb U	24 vials of STA-Desorb U	24 x 15 mL
00973	STA-Cleaner Solution	6 vials of STA-Cleaner Solution	6 x 2.5 L
Disposables			
39430	STA-Satellite Cuvettes	6 rolls of 220 cuvettes and 6 disposable cuvette bins	6 x 220
00741	STA - Microtainer Plastic cups for samples	1 box of 500 Microcontainer	1 x 500
00802	STA-Microcups Microtubes in siliconed glass for reagents control and calibrator plasmas	1 box of 100 STA-Microcups	1 x 100
00797	STA-Mini Reducer Plastic chimneys for reagents (from 4 to 6 mL vials) - Improved stability	1 box of 100 STA-Mini Reducer	1 x 100
00801	STA-Maxi Reducer Plastic chimneys for reagents (from 8 to 15 mL vials) - Improved stability	1 box of 100 STA-Maxi Reducer	1 x 100
27425	White stirring magnet (for STA-Néoplastine and STA-NeoPTimal)	1 item	1
89174	White stirring magnet (for STA-Néoplastine and STA-NeoPTimal)	10 items	1 x 10
26649	Thermal paper	1 roll (25 m)	1
Disposables spare parts			
39356	Equipped needle	1 item	1
26681	5 x 20 T 6.3A fuses	10 fuses	1 x 10
80132	Air filter	3 items	1 x 3
27458	Liquid filter	1 item	1
39738	Liquid filter	10 items	1 x 10

STA SATELLITE MAX / STA SATELLITE USB Auxiliary materials

Cat. NR.	Product name	Product description	Pack.
26555	Ball extractor	1 item	1
88828	Reduction ring DIN 14 (18.5 mm)	2 items	1 x 2
39968	Microcups reducer	2 items	1 x 2
80057	Reducer for microvolume sample	2 items	1 x 2
80091	Reagents carousel n°1	1 item	1
80094	Sample carousel n°1	1 item	1
80092	Reagents carousel n°2	1 item	1
80095	Sample carousel n°2	1 item	1

*Availability depending upon the country.



START MAX

Semi-automated coagulation analyzer, the STart Max offers a user interface standardized with other Max Generation analyzers.

This system manages 4 measuring channels for clotting tests.

Supplied with

- 1 accessory pack

To order*


59017 STart Max (115/230 V)

Parameters*

- PT (Quick and Owren PT)
- APTT
- Fibrinogen
- Thrombin Time
- Reptilase Time
- Extrinsic pathway factors
- Intrinsic pathway factors
- Anti-Xa (coming soon)
- Protein C
- Protein S
- Lupus Anticoagulants
- Calibrators
- Quality controls

Elise M. Project Manager

"Everything has been done to add the modern touches needed, while conserving the simplicity and robustness for which STart has been known for decades."



*Availability depending upon the country.

START MAX / START 4 Disposables

Cat. NR.	Product name	Product description	Pack.
STart Max Disposables			
38876	Cuvettes	150 strips of 4 cuvettes	1
26441	Ball vial	1 vial of 1850 balls	1
86553	Combitips 2.5 mL (position 1 = 50 µL)	1 box of 100	1 x 100
STart 4 Disposables			
38876	Cuvettes	150 strips of 4 cuvettes	1
26441	Ball vial	1 vial of 1850 balls	1
26649	Thermal paper	1 roll (25 m)	1
38051	Finntips 2.5 mL (position 1 = 50 µL)	1 box of 100 Finntips	1 x 100
38052	Finntips 1.25 mL (position 1 = 25 µL)	1 box of 100 Finntips	1 x 100
STart Max & STart 4 Disposables spare parts			
27581	6.3 x 32 T 3A fuses	10 fuses	1 x 10
39478	5 x 20 T 3.15A fuses	10 fuses	1 x 10

START MAX / START 4 Auxiliary materials

Cat. NR.	Product name	Product description	Pack.
STart Max			
86550	Connected pipette	1 item	1
86561	Application notebook STart Max	1 item	1
86562	Notebook STart Max support	1 item	1
86551	Pipette thermal conductor black & grey	2 items (1 black and 1 grey)	1 x 2
26555	Ball extractor	1 item	1
86560	Steel ball dispenser	1 item	1
26605	Reduction ring DIN 14 (18.5 mm)	1 item	1
26610	Reduction ring DIN 18 (22.8 mm)	1 item	1
26729	Reducer for 5 mL vial	1 item	1
26405	Triangular stirring bar (6x12)	1 item	1
89241	10 Triangular stirring bar (6x12)	10 items	1 x 10
89334	USB bar code handheld	1 item	1
86563	Mini USB Adaptor	1 item	1
STart 4			
80109	Connected pipette	1 item	1
26555	Ball extractor	1 item	1
39476	Ball dispenser	1 item	1
26729	Reducer ring for 5 mL tube	1 item	1
26405	Triangular stirring bar (6x12)	1 item	1
26605	Reduction ring DIN 14 (18.5 mm)	1 item	1
26610	Reduction ring DIN 18 (22.8 mm)	1 item	1
38094	Pipette thermal conductor grey	1 item	1
38098	Pipette thermal conductor blue	1 item	1

*Availability depending upon the country.



ST GENESIA SYSTEM

The ST Genesis offers a fully automated system to measure thrombin generation in plasma, with features like the once daily calibration, the reference plasma and the temperature control.

Supplied with

- 1 accessory pack

To order*

86814 ST Genesis

86938 ST Genesis - Research Use Only

Parameters*

- Lag Time
- Peak Height
- Time to Peak
- Velocity Index
- Endogenous Thrombin Potential
- Endogenous Thrombin Potential Inhibition
- Start Tail

Audrey C. Scientific Marketing Manager

"Bringing thrombin generation one step further: automation of the assay was a must-have for the clinical laboratory, and it is now live and accredited in several centers! With many published data and ST Genesis users supporting the idea that thrombin generation is a real aid in diagnosis and patient management, I trust ST Genesis can make TG become a more regular test from the hemostasis lab test panel."




*Availability depending upon the country.

ST GENESIA

Solution & Disposables*

Cat. NR.	Product name	Product description	Pack.
Solution			
00973	STA-Cleaner Solution	6 vials of STA-Cleaner Solution	6 x 2.5 L
Disposables			
86800	STG-Cuvettes	12 racks (40 cuvettes each)	12 x 40
86801	STG-CuvetteWaste	6 items (80 cuvettes capacity each)	1 x 6
86855	Rack for secondary sample tubes ST Genesis	1 item	1
86854	Rack for primary sample tubes ST Genesis	1 item	1
Disposables spare parts			
86853	Kit of 3 filters ST Genesis	3 items	1 x 3
27458	Liquid filter	1 item	1
39738	Liquid filter	10 items	10
39356	Equipped needle	1 item	1
86836	1A ST Genesis fuse	2 items	1 x 2
86837	2A ST Genesis fuse	2 items	1 x 2
86838	3A ST Genesis fuse	2 items	1 x 2
86839	4A ST Genesis fuse	2 items	1 x 2
Disposables - Research Use Only			
86939	STG-Cuvettes	12 racks (40 cuvettes each)	12 x 40

*Availability depending upon the country.



CALIBRATED AUTOMATED THROMBOGRAM

The Calibrated Automated Thrombogram (CAT) is a comprehensive system to measure Thrombin Generation based on fluorescence according to the Hemker & al. method. (Research Product)

Supplied with

- 1 CAT software dedicated to Thrombin Generation measurement
- 1 excitation filter at 390 nM
- 1 emission filter at 460 nM
- 1 user guide
- 1 USB flash
- 1 PC

To order*

86189 CAT System QWERTY

Parameters*

- Lag Time
- Peak Height
- Time to Peak
- Velocity Index
- Endogenous Thrombin Potential
- Start Tail

CAT Disposables

Cat. NR.	Product name	Product description	Pack.
Disposable			
86175	Immulon 2HB Plate	1 box of 50 microtiter plates	1 x 50
Disposables spare parts			
86177	Halogen Lamp	1 item	1
86185	Complete dispensing tube assem.	1 item	1
86182	Reagent tubing	1 item	1
86176	Dispensing tip	1 item	1

*Availability depending upon the country.



PLATELET AGGREGOMETERS TA-8V & TA-4V

The Platelet Aggregometers TA-8V & TA-4V are semi-automated analyzers dedicated to the exploration and evaluation of platelet function in a citrated platelet rich plasma.

Supplied with

- 1 accessory pack
- 1 screen

To order*

58832 TA-8V

58833 TA-4V

PLATELET AGGREGOMETERS Disposable

Cat. NR.	Product name	Product description	Pack.
86921	Glass tubes + stirring bars	1 pack of 1000 pieces	1 x 1000

*Availability depending upon the country.

DIGITAL SOLUTIONS

DIGITAL SOLUTIONS



This catalogue contains information on products which are targeted to a wide range of audiences and could contain product details or information concerning goods which are not accessible or valid in your country.

Coag.One

By Stago

COAG.ONE BY STAGO

Coag.One by Stago is a centralized middleware system that transforms hemostasis laboratory management, and delivers operational efficiency through its ability to organize data, streamline workflow and optimize production from multiple sites.

To order*

Please consult Stago or one of its official distributors.

86765 Coag.One server license

86775 SCE upgrade to Coag.One server license

86780 Coag.One server license Workcell

86795 Coag.One site license

86785 Stago Instrument connection license

89465 Connect.One Gateway kit

86790 SCE upgrade to Coag.One server license Workcell

86770 Coag.One pre-prod server license





ConnectOne

By Stago

CONNECT.ONE BY STAGO

Connect.One offers a secured network architecture allowing remote access and remote monitoring services. Using a single secure entry point, Connect.One is connected devices and send data (quality control results, cost per reportable results, instrument health...) and alarms to Stago managed servers.

To order*

Please consult Stago or one of its official distributors.







MY EXPERT QC

My Expert QC, the real-time peer comparison program by Stago, is a quality control program that continuously monitors the performance of your analyzer based on the routine use of Stago internal quality control (IQC) materials. My Expert QC offers flexible comparison levels from local to global, covering all Stago QC from routine to specialty. My Expert QC program is available 24/7 on www.myexpertqc.com

To order*

Please consult Stago or one of its official distributors.

82127 My Expert QC 1-year contract



KEY FEATURES



- GLOBAL PEER GROUP COMPARISON DEDICATED TO HAEMOSTASIS**
- All routine and specialized tests available
 - Extended WW and nationwide user base
 - Multiple lots of QC available



- REAL-TIME PEER GROUP COMPARISON**
- Automated data transmission
 - Cumulative, monthly or weekly statistical reports
 - Statistical analysis and alerts in seconds



- QC MANAGEMENT MADE EASY**
- Meaningful statistical indicators and graphs
 - Variety of reports and peer groups offered
 - Customizable interface and levels of comparison

IQC
Internal
Quality
Control

Max Generation Systems

QC results are generated
on the analyser

How
does it
work?

STA Coag Expert or Coag.One

Automatic upload of
QC results from Workstation

My Expert QC

Real-time peer group
comparison of IQC results

ConnectOne By Stago

QC results transferred
through secured funnel



MY QUALIRIS QC

Combined with your internal quality control, an External Quality Assessment (EQA) program will allow a comparison of performances with other participants.

Blind plasma covering several concentration levels and reflecting most clinical settings are provided for a regular testing. Peer group comparisons with color-coded alerts and graphs are provided few days after the data submission deadline in monthly reports. Annual reports allow month to month comparison of the results.

My Qualiris QC can be used by all laboratories, regardless of reagents or instrumentation used.

To order*

Please consult Stago
or one of its official distributors.

- 01044 Qualiris QC Premium S1
- 01045 Qualiris QC Premium S2
- 01046 Qualiris QC Experience
- 01047 Qualiris QC Heparin HBPM/LMWH
- 01048 Qualiris QC Heparin HNF/UFH
- 01049 Qualiris QC D-Dimer
- 01050 Qualiris Diagnostic Challenge
- 01063 Qualiris QC Lupus Anticoagulant
- 01187 Qualiris QC Dabigatran
- 01188 Qualiris QC FM-FDP
- 01272 Qualiris QC DOAC Anti-Xa
- 85093 Qualiris QC Premium Factors
- 85094 Qualiris QC Premium Thrombophilia

Pascal M. Qualiris & EIQC Manager

"My Qualiris QC, a constantly evolving independent EQA organization. Stago is adapting to recent developments by offering new FM and DOAC programs. A team dedicated solely to managing Qualiris is always on hand to support you in your daily quality procedures."





REAGENTS



This catalogue contains information on products which are targeted to a wide range of audiences and could contain product details or information concerning goods which are not accessible or valid in your country.

SCREENING TESTS

Extrinsic pathway Prothrombin Time (PT)

Cat. NR.	Product name	Product description	Pack.
Automated reagents			
01163	STA-NeoPTimal 5 Prothrombin Time (PT) (ISI ~ 1)	6 vials of STA-NeoPTimal 5 • 6 vials of solvent	6 x 5 mL
01164	STA-NeoPTimal 10 Prothrombin Time (PT) (ISI ~ 1)	12 vials of STA-NeoPTimal 10 • 12 vials of solvent	12 x 10 mL
01165	STA-NeoPTimal 20 Prothrombin Time (PT) (ISI ~ 1)	12 vials of STA-NeoPTimal 20 • 12 vials of solvent	12 x 20 mL
00665	STA-Néoplastine R 15 Recombinant thromboplastin for determination of Prothrombin Time (PT) (ISI ~ 1)	12 vials of STA-Néoplastine R 15 • 12 vials of solvent	12 x 15 mL
00330	STA-SPA+ Determination of the combined Factors II VII X	12 vials of STA-SPA+ • 12 vials of diluent	12 x 10 mL
Multipurpose reagents			
00374	Néoplastine CI Plus 2 Prothrombin Time (PT) (ISI ~ 1.3)	6 vials of Néoplastine CI Plus 2 • 6 vials of solvent	6 x 2 mL
00122*	SPA 20 Determination of the combined Factors II VII X	12 vials of SPA 20 • 12 vials of diluent	12 x 4 mL

Intrinsic pathway Activated Partial Thromboplastin Time (APTT)

Cat. NR.	Product name	Product description	Pack.
Automated reagents			
00595	STA-PTT Automate 5 Determination of Activated Partial Thromboplastin Time (APTT) (silica activator)	12 vials of STA-PTT Automate 5	12 x 5 mL
00597	STA-C.K. Prest 5 Determination of Activated Partial Thromboplastin Time (APTT) (Kaolin activator)	6 vials of STA-C.K. Prest 5 • 6 vials of solvent + activator	6 x 5 mL
00308	STA-Cephascreen 4 Determination of Activated Partial Thromboplastin Time (APTT) <i>Liquid reagent</i>	12 vials of STA-Cephascreen 4	12 x 4 mL
00310	STA-Cephascreen 10 Determination of Activated Partial Thromboplastin Time (APTT) <i>Liquid reagent</i>	12 vials of STA-Cephascreen 10	12 x 10 mL
Multipurpose reagents			
00598	C.K. Prest 2 Determination of Activated Partial Thromboplastin Time (APTT) (Kaolin activator)	6 vials of C.K. Prest 2 • 6 vials of solvent + activator	6 x 2 mL
00847	C.K. Prest 5 Determination of Activated Partial Thromboplastin Time (APTT) (Kaolin activator)	6 vials of C.K. Prest 5 • 6 vials of solvent + activator	6 x 5 mL
00480	PTT Automate 5 Determination of Activated Partial Thromboplastin Time (APTT) (silica activator)	12 vials of PTT Automate 5	12 x 5 mL
00599	PTT-LA Lupus Anticoagulant APTT based reagent	6 vials of PTT-LA	6 x 2 mL
Reference plasma			
00539	Pool Norm Normal human plasma pool	12 vials of Pool Norm	12 x 1 mL

*For Nordic countries only

*Availability depending upon the country.

Common pathway Fibrinogen

Cat. NR.	Product name	Product description	Pack.
Automated reagent			
00673	STA-Liquid Fib Quantitative determination of fibrinogen by Clauss method <i>Liquid reagent</i>	12 vials of STA-Liquid Fib	12 x 4 mL

Fibrin formation Thrombin Time (TT) - Reptilase Time

Cat. NR.	Product name	Product description	Pack.
Automated & multipurpose reagents			
00611	STA-Thrombin 2 Determination of Thrombin Time	12 vials of STA-Thrombin 2	12 x 2 mL
00669	STA-Thrombin 10 Determination of Thrombin Time	12 vials of STA-Thrombin 10	12 x 10 mL
00614	STA-Reptilase Determination of Reptilase time	6 vials of STA-Reptilase	6 x 2 mL

Magalie L. - Global Solution Product Manager

"1.5 tons of paper saved per year is to contribute to an environmentally responsible action with no compromise on quality of service. A new freely-accessible web platform now ensures the document traceability, where you can download records in all languages and thus reduce the number of documents in paper format."



*Availability depending upon the country.


Combination Reagents

Calibration plasmas-Quality controls

Cat. NR.	Reagent	Calibration plasma	Quality controls
PT			
01163	STA-NeoPTimal 5	Pre-calibrated	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
01164	STA-NeoPTimal 10	Pre-calibrated	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
01165	STA-NeoPTimal 20	Pre-calibrated	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
00665	STA-Néoplastine R 15	Pre-calibrated	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
00330	STA-SPA+	N/A	STA-Routine QC 2 mL - Cat. Nr. 00554
00374	Néoplastine CI Plus 2	Unicalibrator - Cat. Nr. 00625 Etaloquick - Cat. Nr. 00496	Coag Control N+P - Cat. Nr. 00621 System Control N+P - Cat. Nr. 00617
00122*	SPA 20	N/A	STA-Routine QC 2 mL - Cat. Nr. 00554

*For Nordic countries only

*Availability depending upon the country.



Combination Reagents

Calibration plasmas-Quality controls



Cat. NR.	Reagent	Calibration plasma	Quality controls
APTT			
00595	STA-PTT Automate 5	N/A	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
00597	STA-C.K. Prest 5	N/A	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
00308	STA-Cephascreen 4	N/A	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
00310	STA-Cephascreen 10	N/A	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
00598	C.K. Prest 2	N/A	Coag Control N+P - Cat. Nr. 00621 System Control N+P - Cat. Nr. 00617
00847	C.K. Prest 5	N/A	Coag Control N+P - Cat. Nr. 00621 System Control N+P - Cat. Nr. 00617
00480	PTT Automate 5	N/A	Coag Control N+P - Cat. Nr. 00621 System Control N+P - Cat. Nr. 00617
00539	Pool Norm	N/A	STA-Coag Control N+P - Cat. Nr. 00621 System Control N+P - Cat. Nr. 00617
Fibrinogen			
00673	STA-Liquid Fib	Pre-calibrated STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
Thrombin - Reptilase			
00611	STA-Thrombin 2	N/A	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
00669	STA-Thrombin 10	N/A	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554 STA-Routine QC P Plus - Cat. Nr. 00714
00614	STA-Reptilase	N/A	STA-System Control N+P - Cat. Nr. 00678

*Availability depending upon the country.

FACTOR ASSAYS

Deficient plasmas Extrinsic pathway

Cat. NR.	Product name	Product description	Pack.
Automated reagents			
00740	STA-ImmunoDef II Immuno-depleted plasma for factor II assay	6 vials of STA-ImmunoDef II	6 x 1 mL
00745	STA-Deficient II Deficient plasma for factor II assay	6 vials of STA-Deficient II	6 x 1 mL
00744	STA-Deficient V Deficient plasma for factor V assay	6 vials of STA-Deficient V	6 x 1 mL
00743	STA-Deficient VII Deficient plasma for factor VII assay	6 vials of STA-Deficient VII	6 x 1 mL
00738	STA-Deficient X Deficient plasma for factor X assay	6 vials of STA-Deficient X	6 x 1 mL
01113	K-Assay Factor XIII Quantitative determination for factor XIII <i>Liquid reagent</i>	2 vials of buffer reagent • 1 vial of latex	2 x 9.5 mL 1 x 6 mL
01114	K-Assay Factor XIII Calibrator Calibration plasma for factor XIII	5 vials of factor XIII calibrator • 1 vial of diluent	5 x 1 mL 1 x 40 mL
01115	K-Assay Coagulation Control Assayed normal and abnormal control plasmas for factor XIII	5 vials of coagulation control Level 1 • 5 vials of coagulation control Level 2	5 x 2 x 0.5 mL

Sylvie L. - R&D Manager

"The standardization of practices to identify anti-factor antibodies is a challenge for hemostasis laboratories. The use of a complex algorithm, developed in accordance with recommendations by experts in the field, combined with a new graphical representation generally referred to as factor parallelism, will facilitate the interpretation of factor assays."



*Availability depending upon the country.

Deficient plasmas

Intrinsic pathway

Cat. NR.	Product name	Product description	Pack.
Automated reagents			
00728	STA-ImmunoDef VIII Immuno-depleted plasma for factor VIII assay	6 vials of STA-ImmunoDef VIII	6 x 1 mL
T2608*	TriniCHROM Factor VIII:C Quantitative determination of factor VIII:C in human plasma and factor VIII concentrate by chromogenic assay	3 vials of reagent F.IXa • 3 vials of reagent F.X • 3 vials of substrate • 3 vials of dilution buffer	3 x 2 mL 3 x 2 mL 3 x 6 mL 3 x 5 mL
00734	STA-ImmunoDef IX Immuno-depleted plasma for factor IX assay	6 vials of STA-ImmunoDef IX	6 x 1 mL
00759	STA-ImmunoDef XI Immuno-depleted plasma for factor XI assay	6 vials of STA-ImmunoDef XI	6 x 1 mL
00315	STA-ImmunoDef XII Immuno-depleted plasma for factor XII assay	6 vials of STA-ImmunoDef XII	6 x 1 mL
00725	STA-Deficient VIII Immuno-depleted plasma for factor VIII assay	6 vials of STA-Deficient VIII	6 x 1 mL
00724	STA-Deficient IX Immuno-depleted plasma for factor IX assay	6 vials of STA-Deficient IX	6 x 1 mL

Claire M. - R&D Team Leader

"In R&D, we are focused on 3 main activities: developing new products as a first step, but also improving existing reagents and adapting various kits on our instruments. We implement our expertise and know-how on each project in order to offer performant solutions to our customers: for example, automation of TC VIII: C from Tcoag, which permits us to answer to more and more requests for a chromogenic factor VIII determination."



*Availability depending upon the country.

Other methods

Asserachrom - Staclot - TriniCHROM

Cat. NR.	Product name	Product description	Pack.
Factor VII			
00281	Staclot VIIa-rTF Chronometric determination of activated factor VII <i>Research Product</i>	2 vials of deficient plasma VII • 2 vials of buffer • 2 vials of rsTF-Phospholipides • 2 vials of F.VIIa calibrator • 2 vials of control 1 • 2 vials of control 2	2 x 1 mL
00491	Asserachrom VIIa-AT Quantitative determination of factor VIIa-Antithrombin complex by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of sample diluent • 3 vials of anti-AT-peroxidase • 3 vials of TMB • 3 vials of F.VIIa-AT calibrator • 3 vials of F.VIIa-AT control • 3 vials of AT-peroxidase buffer • 1 vial of washing solution	3 x 32 tests
00955	Asserachrom VII:Ag Quantitative determination of factor VII:Ag by ELISA Method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-VII:Ag peroxidase • 3 vials of dilution buffer • 3 vials of F.VII calibrator • 3 vials of washing solution • 3 vials of F.VII control	3 x 32 tests
Factor VIII			
T2608*	TriniCHROM Factor VIII:C Quantitative determination of factor VIII:C in human plasma and factor VIII concentrate by chromogenic assay	3 vials of reagent F.IXa • 3 vials of reagent F.X • 3 vials of substrate • 3 vials of dilution buffer	3 x 2 mL 3 x 2 mL 3 x 6 mL 3 x 5 mL
00280	Asserachrom VIII:Ag Quantitative determination of factor VIII:Ag by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-VIII:Ag peroxidase • 3 vials of dilution buffer • 3 vials of F.VIII:Ag calibrator • 1 vial of washing solution • 3 vials of F.VIII:Ag control	3 x 32 tests
00919	Asserachrom VWF:FVIII B Quantitative determination of the capacity of von Willebrand Factor to bind to factor VIII by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-F.VIII-peroxidase • 3 vials of sample diluent • 3 vials of VWF:F.VIII calibrator • 1 vial of washing solution • 3 vials of VWF:F.VIII control • 3 vials of recomb. F.VIII • 3 vials of anti-F.VIII-peroxidase buffer	3 x 32 tests
Factor IX			
00943	Asserachrom IX:Ag Quantitative determination of factor IX:Ag by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-IX:Ag-peroxidase • 1 vial of dilution buffer • 3 vials of F.IX:Ag calibrator • 1 vial of washing solution • 3 vials of F.IX:Ag control	3 x 32 tests
Factor X			
00956	Asserachrom X:Ag Quantitative determination of factor X:Ag by ELISA Method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-X:Ag-peroxidase • 3 vials of dilution buffer • 3 vials of F.X calibrator • 1 vial of washing solution • 3 vials of F.X control	3 x 32 tests
Factor XIII			
00870	Factor XIII Reagents Factor XIII activity assay	1 vial of purified fibrinogen • 1 vial of kaolin suspension • 1 vial of monochloroacetic acid • 1 vial of calcium thrombin	1 x 2 mL
01113	K-Assay Factor XIII Quantitative determination for factor XIII <i>Liquid reagent</i>	2 vials of buffer reagent • 1 vial of latex	2 x 9.5 mL 1 x 6 mL
01114	K-Assay Factor XIII Calibrator Calibration plasma for factor XIII	5 vials of factor XIII calibrator • 1 vial of diluent	5 x 1 mL 1 x 40 mL
01115	K-Assay Coagulation Control Assayed normal and abnormal control plasmas for factor XIII	5 vials of coagulation control Level 1 • 5 vials of coagulation control Level 2	5 x 2 x 0.5 mL

*Availability depending upon the country.

Combination Reagents

Calibration plasmas-Quality controls



Cat. NR.	Reagent	Calibration plasma	Quality controls
Deficient plasmas			
00740	STA-ImmunoDef II	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00745	STA-Deficient II	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00744	STA-Deficient V	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00743	STA-Deficient VII	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00738	STA-Deficient X	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00728	STA-ImmunoDef VIII	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00734	STA-ImmunoDef IX	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00759	STA-ImmunoDef XI	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00315	STA-ImmunoDef XII	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00725	STA-Deficient VIII	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00724	STA-Deficient IX	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
T2608*	TriniCHROM Factor VIII:C	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678

*Availability depending upon the country.

FIBRIN FORMATION & DEGRADATION

Fibrin formation

Fibrin monomers

Cat. NR.	Product name	Product description	Pack.
Automated reagent			
00543	STA-Liatest FM Quantitative determination of fibrin monomers by immuno turbidimetric method <i>Liquid reagent</i>	6 vials of latex • 6 vials of buffer	6 x 4 mL 6 x 2 mL
Manual reagents			
00857	F.S. Test Detection of soluble fibrin monomers complexes by hemagglutination	4 vials of F.S. Test reagent • 4 vials of positive control • 4 vials of negative control	4 x 0.5 mL
00887	F.S. Test Unit Detection of soluble fibrin monomers complexes by hemagglutination	8 vials of F.S. Test Unit reagent • 8 vials of positive control • 8 vials of negative control • 10 test cards	8 x 0.2 mL
00548	Test cards for F.S. Test and FDP Plasma kits	10 test cards	1 x 10


Degradation products

D-Dimer

Cat. NR.	Product name	Product description	Pack.
Automated reagent			
00662	STA-Liatest D-Di Plus Quantitative determination of D-Dimer levels by immuno-turbidimetric method <i>Liquid reagent</i>	6 vials of latex • 6 vials of buffer	6 x 6 mL 6 x 5 mL
Manual reagents			
00947	Asserachrom D-Di Quantitative determination of D-Dimer levels by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-D Peroxidase • 3 vials of dilution buffer • 3 vials of D-Di calibrator • 1 vial of washing solution • 3 vials of D-Di control	3 x 32 tests
00454	D-Di Test Qualitative and semi quantitative determination of D-Dimer levels by latex agglutination	1 vial of latex • 1 vial of buffer • 1 vial of negative control • 10 test cards • 1 vial of positive control • mixing rods	1 x 1.3 mL
00550	Test cards for D-Di Test kit	10 test cards	1 x 10

Céline D. - Clinical Research Associate

"The DiET study confirmed the excellent clinical performance of Stago D-Di, in accordance with FDA requirements. This is the first prospective clinical study of this scale conducted by Stago (2000 patients, 18 sites ...). Obtaining recognition of the FDA and the publication of the first results in Blood Coagulation and Fibrinolysis (Vol 27, July 2016) is a reward for the whole team."



Degradation products

FDP

Cat. NR.	Product name	Product description	Pack.
Automated reagent			
00649	STA-Liatest FDP Quantitative determination of fibrin and fibrinogen degradation products by immuno-turbidimetric method <i>Liquid reagent</i>	6 vials of latex • 6 vials of buffer	6 x 5 mL 6 x 5 mL
Manual reagents			
00540	FDP Plasma Qualitative and semi quantitative determination of fibrin & fibrinogen degradation products (FDP) in plasma by latex agglutination	1 vial of latex • 1 vial of buffer • 1 vial of negative control • 10 test cards • 1 vial of positive control • mixing rods	1 x 1.3 mL
00541	FDP Plasma (latex)	12 vials of latex	12 x 1.3 mL
00551	FDP Plasma (buffer)	12 vials of buffer	12 x 20 mL
00548	Test cards for F.S. Test and FDP Plasma kits	10 test cards	1 x 10



Combination Reagents

Calibration plasmas-Quality controls

Cat. NR.	Reagent	Calibration plasma	Quality controls
Fibrin formation			
00543	STA-Liatest FM	STA-FM Calibrator - Cat. Nr. 00544	STA-FM Control - Cat. Nr. 00545
Degradation products			
00662	STA-Liatest D-Di Plus	Pre-calibrated	STA-D-Di Control - Cat. Nr. 00868 STA-Liatest Control N+P - Cat. Nr. 00526
00649	STA-Liatest FDP	STA-FDP Calibrator - Cat. Nr. 00655	STA-FDP Control - Cat. Nr. 00654

*Availability depending upon the country.

FIBRINOLYSIS

Fibrinolysis

TAFI - Plasminogen - Antiplasmin - tPA - PAI

Cat. NR.	Product name	Product description	Pack.
00346	STA-Stachrom TAFI Chromogenic assay for the quantitative determination of the Thrombin Activatable Fibrinolysis Inhibitor (TAFI) activity <i>Research Product</i>	2 vials of TAFI activator • 2 vials of substrate • 4 vials of carboxypeptidase A • 2 vials of TAFI calibrator • 2 vials of TAFI control	80 tests
00616*	Asserachrom TAFIa/TAFIai Quantitative determination of activated and/or inactivated TAFI by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-TAFIa/TAFIai Peroxidase • 3 vials of dilution buffer • 3 vials of TAFIa/TAFIai calibrator • 1 vial of washing solution • 3 vials of TAFIa/TAFIai control	3 x 32 tests
00658	STA-Stachrom Plasminogen Chromogenic assay of plasminogen	6 vials of streptokinase • 6 vials of substrate	6 x 3 mL
00659	STA-Stachrom Antiplasmin Chromogenic assay of antiplasmin	4 vials of plasmin • 4 vials of solvent • 4 vials of substrate	4 x 2 mL 4 x 6 mL
00948	Asserachrom tPA Quantitative determination of tissue Plasminogen Activator (tPA) by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-tPA peroxidase • 3 vials of dilution buffer • 3 vials of tPA calibrator • 1 vial of washing solution • 3 vials of tPA control	3 x 32 tests
00807	Stachrom PAI Chromogenic assay of Plasminogen Activator Inhibitor 1 (PAI 1)	2 vials of urokinase • 2 vials of PAI calibrator 1 • 2 vials of plasminogen • 2 vials of PAI calibrator 2 • 2 vials of substrate • 2 vials of PAI calibrator 3	2 x 2 mL
00949	Asserachrom PAI 1 Quantitative determination of Plasminogen Activator Inhibitor 1 (PAI 1) by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-PAI 1 peroxidase • 3 vials of dilution buffer • 3 vials of PAI 1 calibrator • 1 vial of washing solution • 3 vials of PAI 1 control	3 x 32 tests

Combination Reagents

Calibration plasmas-Quality controls



Cat. NR.	Reagent	Calibration plasma	Quality controls
00658	STA-Stachrom Plasminogen	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00659	STA-Stachrom Antiplasmin	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678

*Availability depending upon the country.

ANTICOAGULANT TREATMENTS


Heparins & Direct Anti-Xa

UFH - LMWH - fondaparinux - rivaroxaban - apixaban - edoxaban

Cat. NR.	Product name	Product description	Pack.
Anti-Xa chromogenic assays			
00311	STA-Liquid Anti-Xa 4 Chromogenic assay of UFH, LMWH, fondaparinux (Arixtra), rivaroxaban (Xarelto), apixaban (Eliquis) and edoxaban (Lixiana/Savaysa) <i>Liquid reagent</i>	6 vials of substrate • 6 vials of factor Xa	6 x 4 mL
00322	STA-Liquid Anti-Xa 8 Chromogenic assay of UFH, LMWH, fondaparinux (Arixtra), rivaroxaban (Xarelto), apixaban (Eliquis) and edoxaban (Lixiana/Savaysa) <i>Liquid reagent</i>	6 vials of substrate • 6 vials of factor Xa	6 x 8 mL
00906	Stachrom Heparin Chromogenic assay of UFH and LMWH	4 vials of antithrombin • 4 vials of factor Xa • 2 vial of buffer • 4 vials of substrate	4 x 4 mL
00691	STA-Staclot Heparin 1 Clotting assay of UFH and LMWH	4 vials of substrate plasma • 4 vials of factor Xa • 4 vials of phospholipids Ca ⁺⁺	4 x 1 mL

Tristan H. - Pharmaceutical Development Director

"After 50 years of VKA predominance, direct oral anticoagulants have upset this market and the care of patients. To be closer to this evolution, Stago develops exclusive cooperations with the pharma companies producing each of these drugs in order to make available diagnostic tools to measure their plasma concentration."



*Availability depending upon the country.

HIT

Heparin Induced Thrombocytopenia

Cat. NR.	Product name	Product description	Pack.
01058	STic Expert HIT 5 Lateral Flow Immunoassay to detect IgG antibodies against heparin-PF 4 complexes	5 tests • 1 vial of buffer	1 x 5 tests
01059	STic Expert HIT 20 Lateral Flow Immunoassay to detect IgG antibodies against heparin-PF 4 complexes	20 tests • 1 vial of buffer	1 x 20 tests
00615	Asserachrom HPIA Determination of anti-heparin PF4 antibodies by ELISA method	3 x 2 coated strips • 3 vials of HPIA negative control • 3 vials of HPIA reference • 3 vials of HPIA positive control • 1 vial of washing solution • 6 vials of TMB • 6 vials of anti-human IgG/A/M peroxidase • 3 vials of dilution buffer	3 x 2 x 8 tests
00624	Asserachrom HPIA-IgG Determination of anti-heparin PF4 antibodies of IgG class by ELISA method	3 x 2 coated strips • 3 vials of HPIA negative control • 3 vials of HPIA reference • 3 vials of HPIA positive control • 1 vial of washing solution • 6 vials of TMB • 6 vials of anti-human IgG peroxidase • 3 vials of dilution buffer	3 x 2 x 8 tests

DTI

Direct Thrombin Inhibitors

Cat. NR.	Product name	Product description	Pack.
00992	STA-ECA II Ecarin chromogenic assay of direct thrombin inhibitors (dabigatran)	2 vials of prothrombin • 2 vials of substrate • 2 vials of ecarin	2 x 25 tests

*Availability depending upon the country.

Combination Reagents

Calibration plasmas-Quality controls



Reagent	STA-Liquid Anti-Xa Cat. Nr. 00311 (4 mL) / Cat. Nr. 00322 (8 mL)				
Molecules	UFH & LMWH	fondaparinux	rivaroxaban	apixaban	edoxaban
Calibration type	Dedicated calibration HNF/ HBPM Hybrid calibration	Dedicated calibration	Dedicated calibration	Dedicated calibration	Dedicated calibration
Calibration plasma	STA-Multi-Hep Calibrator Cat. Nr. 00348	STA-Fondaparinux Calibrator Cat. Nr. 00354	STA-Rivaroxaban Calibrator Cat. Nr. 00704	STA-Apixaban Calibrator Cat. Nr. 01075	STA-Edoxaban Calibrator Cat. Nr. 01073
Quality Controls	STA-Quality HNF/ UFH Cat. Nr. 00381 STA-Quality HBPM/ LMWH Cat. Nr. 00686	STA-Fondaparinux Control Cat. Nr. 00355	STA-Rivaroxaban Control Cat. Nr. 00706	STA-Apixaban Control Cat. Nr. 01074	STA-Edoxaban Control Cat. Nr. 01072

Reagent	Stachrom Heparin Cat. Nr. 00906	
Molecules	UFH	LMWH
Calibration type	Dedicated calibration UFH	Dedicated calibration LMWH
Calibration plasma	STA-Hepanorm H Cat. Nr. 00684	STA-Hepanorm HBPM Cat. Nr. 00681
Quality Controls	STA-Heparin Control (UFH) Cat. Nr. 00683	STA-HBPM/LMWH Control Cat. Nr. 00682

Combination Reagents

Calibration plasmas-Quality controls



Cat. NR.	Reagent	Calibration plasma	Quality controls
00992	STA-ECA II	STA-Dabigatran Calibrator - Cat. Nr. 00993	STA-Dabigatran Control - Cat. Nr. 00994

*Availability depending upon the country.

THROMBOPHILIA

Major Physiological Inhibitors

AT - Protein C - APCR - Protein S - C4b-BP

Cat. NR.	Product name	Product description	Pack.
Antithrombin (AT)			
00596	STA-Stachrom AT III 3 Chromogenic assay of Antithrombin	4 vials of thrombin • 4 vials of substrate • 4 vials of thrombin solvent	4 x 3 mL
00672	STA-Stachrom AT III 6 Chromogenic assay of Antithrombin	4 vials of thrombin • 4 vials of substrate • 4 vials of thrombin solvent	4 x 6 mL
00568	Liatest AT III Quantitative determination of Antithrombin by immuno-turbidimetric method	6 vials of latex • 6 vials of buffer	6 x 1 mL
Protein C			
00747	STA-Staot Protein C 1 Clotting assay of Protein C	3 vials of PC deficient plasma • 3 vials of PC activator	3 x 1 mL
00737*	STA-Staot Protein C 3 Clotting assay of Protein C	6 vials of PC deficient plasma • 6 vials of PC activator	6 x 3 mL
00671	STA-Stachrom Protein C Chromogenic assay of Protein C	6 vials of PC activator • 6 vials of substrate	6 x 6 mL
00944	Asserachrom Protein C Quantitative determination of Protein C by ELISA method	3 x 2 coated strips • 3 vials of Protein C control • 1 vial of washing solution • 3 vials of TMB • 3 vials of anti-Protein C peroxidase • 3 vials of dilution buffer • 3 vials of Protein C calibrator	3 x 32 tests
Activated Protein C Resistance (APCR)			
00721	STA-Staot APC R Clotting test for the detection of APC Resistance	4 vials of FV deficient plasma • 4 vials of Control N • 4 vials of APC • 4 vials of Control P • 4 vials of venom	4 x 2 mL
Protein S			
00746	STA-Staot Protein S Clotting assay of Protein S	2 vials of PS deficient plasma • 2 vials of F.Va • 2 vials of APC	2 x 1 mL
00527	STA-Liatest Free Protein S 2 Quantitative determination of Free Protein S by immuno-turbidimetric method <i>Liquid reagent</i>	3 vials of latex • 3 vials of buffer	3 x 2 mL
00516	STA-Liatest Free Protein S 6 Quantitative determination of Free Protein S by immuno-turbidimetric method <i>Liquid reagent</i>	6 vials of latex • 6 vials of buffer	6 x 6 mL
00945	Asserachrom Total Protein S Quantitative determination of Total Protein S by ELISA method	3 x 2 coated strips • 3 vials of Total Protein S calibrator • 3 vials of dilution buffer • 3 vials of TMB • 3 vials of Total Protein S control • 1 vial of washing solution • 3 vials of anti-Total Protein S peroxidase	3 x 32 tests
00946	Asserachrom Free Protein S Quantitative determination of Free Protein S by ELISA method	3 x 2 coated strips • 3 vials of Free Protein S calibrator • 3 vials of dilution buffer • 3 vials of TMB • 3 vials of Free Protein S control • 1 vial of washing solution • 3 vials of anti-Free Protein S peroxidase	3 x 32 tests
C4b-BP			
00581	Liatest C4b-BP Quantitative determination of C4b-BP by immuno-turbidimetric method Research Product	6 vials of latex • 6 vials of buffer	6 x 1 mL

*Availability depending upon the country.

Antiphospholipid syndrome Lupus Anticoagulants (LA)

Cat. NR.	Product name	Product description	Pack.
DRVV Screening and Confirmation			
00339	STA-Sta clot DRVV Screen 2 Detection of the LA by the diluted Russel's viper venom test	12 vials of STA-Sta clot DRVV Screen 2	12 x 2 mL
00333	STA-Sta clot DRVV Screen 5 Detection of the LA by the diluted Russel's viper venom test	12 vials of STA-Sta clot DRVV Screen 5	12 x 5 mL
00334	STA-Sta clot DRVV Confirm Confirmation of LA by the diluted Russel's viper venom test	12 vials of STA-Sta clot DRVV Confirm	12 x 2 mL
APTT Method Screening and Confirmation			
00599	PTT-LA Lupus Anticoagulant APTT based reagent	6 vials of PTT-LA	6 x 2 mL
00600	Sta clot LA Hexagonal phospholipids screening and confirmatory assay for LA	2 vials of hexagonal phospholipids • 2 vials of buffer • 4 vials of normal plasma • 2 vials of solvent • 2 vials of PTT LS	10 tests
Reference plasma			
00539	Pool Norm Normal human plasma pool	12 vials of Pool Norm	12 x 1 mL

Other Inhibitors HCII - TFPI - sEPCR

Cat. NR.	Product name	Product description	Pack.
Heparin Cofactor II (HCII)			
00851	Stachrom HCII Chromogenic assay of Heparin Cofactor II <i>Research Product</i>	6 vials of thrombin • 6 vials of substrate • 2 vials of buffer	6 x 2 mL
Tissue Factor Pathway Inhibitor (TFPI)			
00261	Asserachrom Total TFPI Quantitative determination of Total Tissue Factor Pathway Inhibitor (TFPI) by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of anti-Total TFPI peroxidase • 1 vial of washing solution • 6 tablets of OPD • 3 vials of Total TFPI calibrator • 6 tablets of urea peroxide • 3 vials of Total TFPI control • 1 vial of dilution buffer	3 x 32 tests
Soluble Endothelial Protein C Receptor (sEPCR)			
00264	Asserachrom sEPCR Quantitative determination of soluble Endothelial Protein C Receptor by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 1 vial of washing solution • 3 vials of dilution buffer • 3 vials of sEPCR calibrator • 3 vials of sEPCR control • 3 vials of anti-sEPCR peroxidase	3 x 32 tests

*Availability depending upon the country.

Combination Reagents

Calibration plasmas-Quality controls



Cat. NR.	Reagent	Calibration plasma	Quality controls
Antithrombin (AT)			
00596	STA-Stachrom AT III 3	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554
00672	STA-Stachrom AT III 6	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678 STA-Coag Control N+P - Cat. Nr. 00679 STA-Routine QC 2 mL - Cat. Nr. 00554
Protein C			
00747	STA-Staclot Protein C 1	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00737*	STA-Staclot Protein C 3	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00671	STA-Stachrom Protein C	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
Protein S			
00746	STA-Staclot Protein S	STA-Unicalibrator - Cat. Nr. 00675	STA-System Control N+P - Cat. Nr. 00678
00527	STA-Liatest Free Protein S 2	Pre-calibrated	STA-Liatest Control N+P - Cat. Nr. 00526
00516	STA-Liatest Free Protein S 6	Pre-calibrated	STA-Liatest Control N+P - Cat. Nr. 00526
LA			
00339	STA-Staclot DRVV Screen 2	N/A	STA-Control LA 1+2 - Cat. Nr. 00201
00333	STA-Staclot DRVV Screen 5	N/A	STA-Control LA 1+2 - Cat. Nr. 00201
00334	STA-Staclot DRVV Confirm	N/A	STA-Control LA 1+2 - Cat. Nr. 00201
00600	Staclot LA	N/A	STA-Control LA 1+2 - Cat. Nr. 00201



*Availability depending upon the country.

PRIMARY HEMOSTASIS

Von Willebrand Disease Von Willebrand Factor (VWF)

Cat. NR.	Product name	Product description	Pack.
VWF			
01191*	STA-VWF:Rco Automated assay for the determination of Ristocetin Cofactor activity of Von Willebrand Factor	3 vials of human lyophilized platelets • 3 vials of Ristocetin • 1 vials of TBS buffer • 2 vials of calibrator • 3 vials of normal control • 3 vials of abnormal control	3 x 20 tests
00518	STA-Liatest VWF:Ag Quantitative determination of VWF by immuno-turbidimetric method	4 vials of latex • 4 vials of buffer • 4 vials of latex diluent	4 x 5 mL
00942	Asserachrom VWF:Ag Quantitative determination of VWF by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-VWF peroxidase • 1 vial of dilution buffer • 3 vials of VWF calibrator • 1 vial of washing solution • 3 vials of VWF control	3 x 32 tests
00239	Asserachrom VWF:CB Quantitative determination of the capacity of VWF to bind to Collagen by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-VWF peroxidase • 3 vials of dilution buffer • 3 vials of VWF:CB calibrator • 1 vial of washing solution • 3 vials of VWF:CB control	3 x 32 tests
00919	Asserachrom VWF:FVIII B Quantitative determination of the capacity of VWF to bind to F.VIII by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-F.VIII peroxidase • 3 vials of dilution buffer • 3 vials of VWF:F.VIII calibrator • 1 vial of washing solution • 3 vials of VWF:F.VIII control • 3 vials of recomb. F.VIII • 3 vials of recomb. F.VIII and peroxidase buffer	3 x 32 tests
00501	Ristocetin <i>Research Product</i>	1 vial of Ristocetin	1 x 5 mL

Combination Reagents Calibration plasmas-Quality controls



Cat. NR.	Reagent	Calibration plasma	Quality controls
VWF Factor			
00518	STA-Liatest VWF:Ag	STA-Liatest VWF:Ag Calibrator - Cat. Nr. 00520	STA-Liatest Control N+P - Cat. Nr. 00526

*Availability depending upon the country.

Platelets

Flow Cytometry & ELISA*

Cat. NR.	Product name	Product description	Pack.
Flow Cytometry assays			
00449	PLT VASP/P2Y12 For the monitoring of P2Y12 ADP receptor antagonists by flow cytometry	1 vial of diluent • 1 vial of PGE1 • 1 vial of PGE1 + ADP • 1 vial of fixative agent • 1 vial of anti-VASP P mouse monoclonal antibody • 1 vial of negative isotypic control • 1 vial of staining reagent	10 samples
00111	PLT GP Receptors Quantitative determination of platelet surface glycoproteins by flow cytometry	1 vial of diluent • 1 vial of fixative agent • 1 vial of negative isotypic control • 1 vial of TRAP • 1 vial of anti-GPIIb/IIIa • 1 vial of anti-GPIb • 1 vial of anti-GPIIIa • 1 vial of anti-GMP 140 • 1 vial of staining reagent • 1 vial of calibrator	5 samples
00452	Platelet PAIg Kit for Platelet Associated Immunoglobulin quantitation by flow cytometry <i>Research Product</i>	1 vial of diluent • 1 vial of calibrator • 1 vial of Mab • 1 vial of staining reagent • 1 vial of negative isotypic control • 1 vial of buffer	10 samples
00457	Platelet Calibrator Kit for customized platelet antigen quantitation by flow cytometry <i>Research Product</i>	1 vial of diluent • 1 vial of negative isotypic control IgG2b • 1 vial of negative isotypic control IgG1 • 1 vial of calibrator • 1 vial of negative isotypic control IgG2a • 1 vial of staining reagent	50 samples
00418	Platelet GP Screen Kit for customized platelet glycoprotein quantitation by flow cytometry <i>Research Product</i>	1 vial of diluent • 1 vial of Mab2 anti-GPIa • 1 vial of calibrator • 1 vial of Mab2 anti-GPIb • 1 vial of staining reagent • 1 vial of Mab2 anti-GPIIIa	10 samples
00420	Megamix Beads for cytometer settings in microparticles analysis <i>Research Product</i>	1 vial of beads	50 tests
01077	Megamix-Plus FSC Beads for FSC-optimized cytometer settings in microparticles analysis <i>Research Product</i>	1 vials of beads	50 tests
01078	Megamix-Plus SSC Beads for SSC-optimized cytometer settings in microparticles analysis <i>Research Product</i>	1 vial of beads	50 tests
01169	MP-Count Beads Absolute counting microparticles by flow cytometry on PPP or purified microparticles <i>Research Product</i>	1 vial of 3 mL	100 tests
ELISA assay			
01076	Cy-Quant VASP/P2Y12 For the monitoring of P2Y12 ADP receptor antagonists by ELISA	96 divisible anti-VASP coated wells • 1 vial of washing solution • 3 vials of PGE1 • 1 vial of dilution buffer • 3 vials of PGE1 + ADP • 1 vial of TMB • 1 vial of anti-VASP-P peroxidase • 1 vial of stop solution • 1 vial of lysis buffer	96 unitary tests
01514	Cy-Quant MV/TF Activity Kit for quantitative measurement of Tissue Factor activity <i>Research Product</i>	1 vial of washing solution • 1 vial of non blocking TF Mab • 1 vial of anti-TF Mab • 2 vials of FVII + FX mix • 1 vial of dilution buffer • 1 vial of stop solution • 1 vial of chromogenic substrate • 2 vials of Calibrator • 2 vials of low activity control • 2 vials of high activity control	96 tests

*Availability depending upon the country.

Platelets

Platelet aggregation and activation markers*

Cat. NR.	Product name	Product description	Pack.
86922	Arachidonic acid	3 vials of Arachidonic acid	3 x 1 mL
86923	ADP	3 vials of ADP	3 x 1 mL
86924	Collagen	3 vials of Collagen	3 x 1 mL
86925	Epinephrin	3 vials of Epinephrin	3 x 1 mL
86926	TRAP 6	3 vials of TRAP 6	3 x 1 mL
00501	Ristocetin <i>Research Product</i>	1 vial of Ristocetin	1 x 5 mL

Activation markers

00950	Asserachrom β TG Quantitative determination of β Thromboglobulin (β TG) by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti- β TG peroxidase • 3 vials of dilution buffer • 3 vials of β TG calibrator • 1 vial of washing solution • 3 vials of β TG control	3 x 32 tests
00951	Asserachrom PF4 Quantitative determination of Platelet Factor 4 (PF4) by ELISA method	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-PF4 peroxidase • 3 vials of dilution buffer • 3 vials of PF4 calibrator • 1 vial of washing solution • 3 vials of PF4 control	3 x 32 tests

*Availability depending upon the country.

Monoclonal Antibodies

Flow Cytometry*

Cat. NR.	Product name	Product description	Pack.
Anti-platelet markers			
01033	CD32, clone 2B2, purif.		0.1 mg
01031	CD36, clone 10.5, purif.		0.1 mg
01030	CD36, clone 10.5, FITC		100 tests
01032	CD36, clone 10.5, PE		100 tests
01025	CD41, clone PL2 49, purif.		0.1 mg
01024	CD41, clone PL2 49, FITC		100 tests
01026	CD41, clone PL2 49, PE		100 tests
01028	CD42b, clone ALMA 19, purif.		0.1 mg
01027	CD42b, clone ALMA 19, FITC		100 tests
01029	CD42b, clone ALMA 19, PE		100 tests
01017	CD61, clone 4F8, purif.		0.1 mg
01016	CD61, clone 4F8, FITC		100 tests
01022	CD62P, clone LYP20, purif.		0.1 mg
01021	CD62P, clone LYP20, FITC		100 tests
01023	CD62P, clone LYP20, PE		100 tests
01005	Fibrinogen, clone 9F9, FITC		100 tests
01083	GPVI, clone 1G5, purif.		0.1 mg
01084	GPVI, clone 1G5, PE		100 tests
Anti-endothelial cell markers			
01006	CD146, clone COM3D9, purif.		100 tests
01007	CD146, clone COM2F6, purif.		100 tests
01008	CD146, clone COM5G6, purif.		100 tests
01148	CD146, clone COM7A4, purif.		0.1 mg
01010	CD146, clone COM7A4, FITC		100 tests
01009	CD146, clone COM7A4, Biot.		100 tests
01149	CD146, clone S-ENDO 1, purif.		0.1 mg
01013	CD146, clone S ENDO 1, FITC		100 tests
01015	CD146, clone S ENDO 1, PE		100 tests
01012	CD146, clone S ENDO 1, Biot.		100 tests

*Availability depending upon the country.

Monoclonal Antibodies

Flow Cytometry*

Cat. NR.	Product name	Product description	Pack.
Isotypic negative controls			
01019	Ctl. neg. IgG1 purif. (2DNP2H11)		0.1 mg
01018	Ctl. neg. IgG1 FITC (2DNP2H11)		100 tests
01020	Ctl. neg. IgG1 PE (2DNP2H11)		100 tests
01038	Ctl. neg. IgG2a purif. (2DNP16C12)		0.1 mg
01037	Ctl. neg. IgG2a FITC (2DNP16C12)		100 tests
01039	Ctl. neg. IgG2a PE (2DNP16C12)		100 tests
01035	Ctl. neg. IgG2b purif. (2DNP14G5)		0.1 mg
01034	Ctl. neg. IgG2b FITC (2DNP14G5)		100 tests
01036	Ctl. neg. IgG2b PE (2DNP14G5)		100 tests



*Availability depending upon the country.

GLOBAL TESTS

Thrombin Generation ST Genesia*

Cat. NR.	Product name	Product description	Pack.
01277	STG-BleedScreen <i>Research Product</i>	3 vials of STG-BleedScreen • 3 vials of STG-RefPlasma BLS • 3 vials of STG-QualiTest Norm BLS • 3 vials of STG-QualiTest Low BLS	3 x 4 x 1 mL
01278	STG-DrugScreen <i>Research Product</i>	3 vials of STG-DrugScreen • 3 vials of STG-RefPlasma DS • 3 vials of STG-QualiTest Norm DS • 3 vials of STG-QualiTest Low DS	3 x 4 x 1 mL
01133	STG-ThromboScreen	3 vials of STG-ThromboScreen -TM • 3 vials of STG-ThromboScreen+TM • 3 vials of STG-RefPlasma TS • 3 vials of STG-QualiTest High TS • 3 vials of STG-QualiTest Norm TS • 3 vials of STG-QualiTest Low TS	3 x 6 x 1 mL
01141	STG-Cal&Fluo	3 vials of STG-ThrombiCal • 3 vials of STG-FluoStart • 3 vials of STG-FluoSet	3 x 2 mL 3 x 1.5 mL 3 x 1.5 mL
01140	STG-ThrombiClean	6 vials of STG-ThrombiClean	6 x 2 mL



*Availability depending upon the country.

CALIBRATORS & CONTROLS

Calibration plasmas

Automated & Multipurpose reagents

Cat. NR.	Product name	Product description	Pack.
Automated reagents			
00675	STA-Unicalibrator Calibration plasma for PT (%), fibrinogen, factor assays (II, V, VII, VIII, IX, X, XI, XII), AT, PC, PS, plasminogen, antiplasmin (activity)	6 vials of STA-Unicalibrator	6 x 1 mL
00520	STA-VWF:Ag Calibrator Calibration plasma for assay of von Willebrand factor by immuno turbidimetric method	6 vials of STA-VWF:Ag Calibrator	6 x 1 mL
00544	STA-FM Calibrator Calibration plasmas for assay of fibrin monomers by immuno turbidimetric method (STA-Liatest FM)	2 vials of STA-FM Calibrator 1 • 2 vials of STA-FM Calibrator 2 • 2 vials of STA-FM Calibrator 3 • 2 vials of STA-FM Calibrator 4 • 2 vials of STA-FM Calibrator 5	2 x 5 x 1 mL
00655	STA-FDP Calibrator Calibration plasmas for assay of fibrin and fibrinogen degradation products by immuno-turbidimetric method (STA-Liatest FDP)	2 vials of STA-FDP Calibrator 1 • 2 vials of STA-FDP Calibrator 2 • 2 vials of STA-FDP Calibrator 3 • 2 vials of STA-FDP Calibrator 4 • 2 vials of STA-FDP Calibrator 5	2 x 5 x 1 mL
00684	STA-Hepanorm H Calibration plasmas for UFH assays using anti-Xa method (STA-Staclot Heparin/Stachrom Heparin)	4 vials of STA-Hepanorm H 0 • 4 vials of STA-Hepanorm H 3 • 4 vials of STA-Hepanorm H 6	4 x 3 x 1 mL
00681	STA-Hepanorm HBPM/LMWH Calibration plasmas for LMWH assays using anti-Xa method (STA-Staclot Heparin/Stachrom Heparin)	4 vials of STA-Hepanorm HBPM/LMWH 0 • 4 vials of STA-Hepanorm HBPM/LMWH 4 • 4 vials of STA-Hepanorm HBPM/LMWH 9	4 x 3 x 1 mL
00348	STA-Multi Hep Calibrator Calibration plasmas for heparin assays (UFH and LMWH) using anti-Xa method (STA-Liquid Anti-Xa)	4 vials of STA-Multi Hep Calibrator 0 • 4 vials of STA-Multi Hep Calibrator 4 • 4 vials of STA-Multi Hep Calibrator 7 • 4 vials of STA-Multi Hep Calibrator 10 • 4 vials of STA-Multi Hep Calibrator 18	4 x 5 x 1 mL
00354	STA-Fondaparinux Calibrator Calibration plasmas for fondaparinux (Arixtra) assay using anti-Xa method (STA-Liquid Anti-Xa)	2 vials of STA-Fondaparinux Calibrator 0 • 2 vials of STA-Fondaparinux Calibrator 9 • 2 vials of STA-Fondaparinux Calibrator 18	2 x 3 x 1 mL
00704	STA-Rivaroxaban Calibrator Calibration plasmas for rivaroxaban (Xarelto) assay using anti-Xa method (STA-Liquid Anti-Xa)	3 vials of STA-Rivaroxaban 0 • 3 vials of STA-Rivaroxaban 1 • 3 vials of STA-Rivaroxaban 2 • 3 vials of STA-Rivaroxaban 3	3 x 4 x 1 mL
01075	STA-Apixaban Calibrator Calibration plasmas for apixaban (Eliquis) assay using anti-Xa method (STA-Liquid Anti-Xa)	3 vials of STA-Apixaban Calibrator 0 • 3 vials of STA-Apixaban Calibrator 1 • 3 vials of STA-Apixaban Calibrator 2 • 3 vials of STA-Apixaban Calibrator 3	3 x 4 x 1 mL
01073	STA-Edoxaban Calibrator Calibration plasmas for edoxaban (Lixiana, Savaysa) assay using anti-Xa method (STA-Liquid Anti-Xa)	3 vials of STA-Edoxaban Calibrator 0 • 3 vials of STA-Edoxaban Calibrator 1 • 3 vials of STA-Edoxaban Calibrator 2 • 3 vials of STA-Edoxaban Calibrator 3	3 x 4 x 1 mL
00993	STA-Dabigatran Calibrator Calibration plasmas for dabigatran (Pradaxa) assay using ecarin chromogenic method (STA-ECA II)	3 vials of STA-Dabigatran Calibrator 0 • 3 vials of STA-Dabigatran Calibrator 1 • 3 vials of STA-Dabigatran Calibrator 2 • 3 vials of STA-Dabigatran Calibrator 3 • 3 vials of STA-Dabigatran Calibrator 4	2 x 5 x 1 mL
Multipurpose reagents			
00625	Unicalibrator Calibration plasma for PT (%), fibrinogen, factor assays (II, V, VII, VIII, IX, X, XI, XII), AT (antigen), PC, PS (activity)	6 vials of Unicalibrator	6 x 1 mL
00496	Etaloquick Calibration plasmas for the determination of PT (%), INR)	4 vials of Etaloquick 1 • 4 vials of Etaloquick 2 • 4 vials of Etaloquick 3	4 x 3 x 0.5 mL

*Availability depending upon the country.

Control plasmas

Automated reagents*

Cat. NR.	Product name	Product description	Pack.
Automated reagents			
00554	STA-Routine QC 2 mL Normal and abnormal control plasmas for PT/INR, Owren's PT, APTT, fibrinogen, TT (N) and AT (activity) (stable 24h)	12 vials of STA-Routine QC 2 mL N • 12 vials of STA-Routine QC 2 mL P	12 x 2 x 2 mL
00714	STA-Routine QC P Plus Abnormal control plasma for PT, APTT, fibrinogen, TT (N) and AT (activity)	24 vials of STA-Routine QC P Plus	24 x 2 mL
00679	STA-Coag Control N+P Normal and abnormal control plasmas for PT, APTT, fibrinogen, TT (N) and AT (activity)	12 vials of STA-Coag Control N • 12 vials of STA-Coag Control P	12 x 2 x 1 mL
00678	STA-System Control N+P Normal and abnormal control plasmas for PT, APTT, fibrinogen, TT (N), reptilase time (N), factor assays (II, V, VII, VIII, IX, X, XI, XII), AT, PC, PS, plasminogen, antiplasmin (activity)	12 vials of STA-System Control N • 12 vials of STA-System Control P	12 x 2 x 1 mL
00677	STA-Coag Control N ABN Plus (2 mL) Normal and pathological control plasmas for PT/INR, Owren's PT, APTT, fibrinogen, TT (N) and AT (activity) (stable 24h)	12 vials of STA-Coag Control N ABN Plus • 12 vials of STA-Coag Control ABN Plus	12 x 2 x 2 mL
00526	STA-Liatest Control N+P Normal and abnormal control plasmas for assays of VWF, free Protein S and D-Dimer by immuno turbidimetric method	12 vials of STA-Liatest Control N • 12 vials of STA-Liatest Control P	12 x 2 x 1 mL
00868*	STA-D-Di Control Assayed liquid controls for D-Dimer tests by immuno-turbidimetric method (STA-Liatest D-Di, STA-Liatest D-Di Plus) <i>Liquid reagent</i>	6 x 2 vials of STA-D-Di Control 1 • 6 x 2 vials of STA-D-Di Control 2	6 x 2 x 2 mL
00654	STA-FDP Control Control plasmas for assay of fibrin and fibrinogen degradation products by immuno turbidimetric method (STA-Liatest FDP)	6 vials of STA-FDP Control 1 • 6 vials of STA-FDP Control 2	6 x 2 x 1.5 mL
00545	STA-FM Control Control plasmas for assay of fibrin monomers by immuno turbidimetric method (STA-Liatest FM)	6 vials of STA-FM Control 1 • 6 vials of STA-FM Control 2	6 x 2 x 1 mL
00201	STA-Control LA 1+2 Control plasmas for LA tests	3 vials of STA-Control LA 1 • 3 vials of STA-Control LA 2	3 x 2 x 1 mL
00682	STA-HBPM/LMWH Control Control plasmas for LMWH assays using anti-Xa method (STA-Staclot Heparin/Stachrom Heparin)	6 vials of STA-Quality HBPM/LMWH Control 3 • 6 vials of STA-Quality HBPM/LMWH Control 8	6 x 2 x 1 mL
00683	STA-Heparin Control Control plasmas for UFH assays using anti-Xa method (STA-Staclot Heparin/Stachrom Heparin)	6 vials of STA-Heparin Control 2 • 6 vials of STA-Heparin Control 5	6 x 2 x 1 mL
00381	STA-Quality HNF/UFH Control plasmas for UFH assay using anti-Xa method (STA-Liquid Anti-Xa)	6 vials of STA-Quality HNF/UFH 2 • 6 vials of STA-Quality HNF/UFH 7	6 x 2 x 1 mL
00686	STA-Quality HBPM/LMWH Control plasmas for LMWH assay using anti-Xa method (STA-Liquid anti-Xa)	6 vials of STA-Quality HBPM/LMWH 8 • 6 vials of STA-Quality HBPM/LMWH 14	6 x 2 x 1 mL
00355	STA-Fondaparinux Control Control plasmas for Fondaparinux (Arixtra) Assay using anti-Xa method (STA-Liquid Anti-Xa)	3 vials of STA-Fondaparinux Control 5 • 3 vials of STA-Fondaparinux Control 14	3 x 2 x 1 mL

*Availability depending upon the country.

Control plasmas

Automated reagents*

Cat. NR.	Product name	Product description	Pack.
00706	STA-Rivaroxaban Control Control plasmas for Rivaroxaban (Xarelto) assay using anti-Xa method (STA-Liquid Anti-Xa)	3 vials of STA-Rivaroxaban Control 1 • 3 vials of STA-Rivaroxaban Control 2	3 x 2 x 1 mL
01074	STA-Apixaban Control Control plasmas for apixaban (Eliquis) assay using anti-Xa method (STA-Liquid Anti-Xa)	3 vials of STA-Apixaban Control 1 • 3 vials of STA-Apixaban Control 2	3 x 2 x 1 mL
01072	STA-Edoxaban Control Control plasmas for edoxaban (Lixiana, Savaysa) assay using anti-Xa method (STA-Liquid Anti-Xa)	3 vials of STA-Edoxaban Control 1 • 3 vials of STA-Edoxaban Control 2	3 x 2 x 1 mL
00994	STA-Dabigatran Control Control plasmas for dabigatran (Pradaxa) assay using ecarin chromogenic method (STA-ECA II)	3 vials of STA-Dabigatran Control 1 • 3 vials of STA-Dabigatran Control 2	3 x 2 x 1 mL

Control plasmas

Multipurpose reagents

Cat. NR.	Product name	Product description	Pack.
00621	Coag Control N+P Normal and abnormal control plasmas for PT, APTT, fibrinogen, TT (N)	12 vials of Coag Control N • 12 vials of Coag Control P	12 x 2 x 1 mL
00617	System Control N+P Normal and abnormal control plasmas for PT, APTT fibrinogen, TT (N), reptilase time (N), factor assays (II, V, VII, VIII, IX, X, XI, XII), AT (antigen), PC/PS (activity)	12 vials of System Control N • 12 vials of System Control P	12 x 2 x 1 mL

STic Expert HIT

Cat. NR.	Product name	Product description	Pack.
01513*	STic Expert HIT Control Positive and negative control plasmas for STic Expert HIT	1 vial of Negative Control • 1 vial of Positive Control	1 x 2 x 0.25mL

Correlation Plasma sets

ExpertCor

Cat. NR.	Product name	Product description	Pack.
01494	ExpertCor Routine <i>Research Product</i>	30 vials for PT, APTT and fibrinogen • 5 vials for INR	35 x 1.5 mL
01307	ExpertCor DDI <i>Research Product</i>	20 vials	20 x 1 mL
01308	ExpertCor UFH <i>Research Product</i>	10 vials	10 x 1 mL
01309	ExpertCor LMWH <i>Research Product</i>	10 vials	10 x 1 mL

*Availability depending upon the country.

External Quality Control Program

My Qualiris QC

Cat. NR.	Product name	Product description	Pack.
01044	Qualiris QC Premium S1 Unassayed plasmas for external quality assessment of PT, APTT, Fibrinogen, TT, Reptilase Time, Endogenous Factors*, Exogenous Factors*, Factor XIII*, VWF*, AT*, PC*, PS*, Antiplasmin*, Plasminogen*	12 vials of plasma - Semester 1	12 x 1 mL
01045	Qualiris QC Premium S2 Unassayed plasmas for external quality assessment of PT, APTT, Fibrinogen, TT, Reptilase Time, Endogenous Factors*, Exogenous Factors*, Factor XIII*, VWF*, AT*, PC*, PS*, Antiplasmin*, Plasminogen*	12 vials of plasma - Semester 2	12 x 1 mL
01046**	Qualiris QC Experience Unassayed plasmas for external quality assessment of PT, APTT, Fibrinogen, TT, Reptilase Time, Endogenous Factors*, Exogenous Factors*, Factor XIII*, VWF*, AT*, PC*, PS*, Antiplasmin*, Plasminogen*	6 vials of plasma - Full year	6 x 1 mL
01049	Qualiris QC D-Dimer Unassayed plasmas for external quality of D-Dimer	6 vials of plasma - Full year	6 x 1 mL
01048	Qualiris QC Heparin HNF/UFH Unassayed plasmas for external quality of UFH	6 vials of plasma - Full year	6 x 1 mL
01047	Qualiris QC Heparin HBPM/LMWH Unassayed plasmas for external quality assessment of LMWH	6 vials of plasma - Full year	6 x 1 mL
01063	Qualiris QC Lupus Anticoagulant Unassayed plasmas for external quality assessment of Lupus Anticoagulant	6 vials of plasma - Full year	6 x 1 mL
01050	Qualiris Diagnostic Challenge Unassayed plasmas with 3 associated clinical case studies	6 vials of plasma - Full year	2 x 3 x 1 mL
01188	Qualiris QC FM-FDP Unassayed plasmas for external quality assessment of FM-FDP	6 vials of plasma - Full year	2 x 3 x 1 mL
01272	Qualiris QC DOAC Anti-Xa Unassayed plasmas for external quality assessment of DOAC Anti-Xa (rivaroxaban, apixaban & edoxaban)	9 vials of plasma - Full year	3 x 3 x 3 mL
01187	Qualiris QC Dabigatran Unassayed plasmas for external quality assessment of Dabigatran	4 vials of plasma - Full year	2 x 2 x 1 mL
85093	Qualiris QC Premium Factors Monthly external quality assessment for factors II, V, VII, VIII, IX, X, XI, XII, XIII, VWF	Factor Option, utilise the Qualiris QC - Premium vials of plasma	N/A
85094	Qualiris QC Premium Thrombophilia Monthly external quality assessment of Antithrombin, Protein C, Protein S, Plasminogen and Antiplasmin	Thrombophilia Option, utilise the Qualiris QC - Premium vials of plasma	N/A

*Option available depending on the program subscription

**Availability depending upon the country.

OTHER REAGENTS

Auxiliary Reagents Solutions & Test Cards

Cat. NR.	Product name	Product description	Pack.
Solutions			
00975	STA-Desorb U	24 vials of STA-Desorb U	24 x 15 mL
00973	STA-Cleaner Solution	6 vials of STA-Cleaner Solution	6 x 2.5 L
00367	STA-CaCl ₂ 0.025 M	24 vials of STA-CaCl ₂ 0.025 M	24 x 15 mL
00360	STA-Owren Koller	24 vials of STA-Owren Koller • Buffer • pH 7.35	24 x 15 mL
00279	PEG 25%	6 vials of PEG 25%	6 x 2.5 mL
00555*	Asserachrom Washing Solution	2 vials of washing solution for Asserachrom	2 x 50 mL
Test Cards			
00550	Test cards for D-Di Test kit	10 test cards	1 x 10
00548	Test cards for F.S. Test and FDP Plasma kits	10 test cards	1 x 10

*Availability depending upon the country.

RESEARCH PRODUCTS



This catalogue contains information on products which are targeted to a wide range of audiences and could contain product details or information concerning goods which are not accessible or valid in your country.

RESEARCH PRODUCTS

Research kits

Cat. NR.	Product name	Product description	Pack.
Activity methods			
00281	Sta clot VIIa-rTF Chronometric determination of activated factor VII <i>Research Product</i>	2 vials of deficient plasma VII • 2 vials of buffer • 2 vials of rsTF-Phospholipides • 2 vials of F.VIIa calibrator • 2 vials of control 1 • 2 vials of control 2	2 x 1 mL
00851	Stachrom HCII Chromogenic assay of Heparin Cofactor II <i>Research Product</i>	6 vials of thrombin • 6 vials of substrate • 2 vials of buffer	6 x 2 mL
00346	STA-Stachrom TAFI Chromogenic assay for the quantitative determination of the Thrombin Activatable Fibrinolysis Inhibitor (TAFI) activity <i>Research Product</i>	2 vials of TAFI activator • 2 vials of substrate • 4 vials of carboxypeptidase A • 2 vials of TAFI calibrator • 2 vials of TAFI control	80 tests
00429	STA-Procoag PPL Chronometric determination of procoagulant phospholipid activity <i>Research Product</i>	3 vials of procoagulant phospholipid depleted plasma • 3 vials of F. Xa • 3 vials of Control N • 3 vials of Control P	3 x 40 tests
ELISA method			
00955	Asserachrom VII:Ag Quantitative determination of factor VII:Ag by ELISA Method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-VII:AG-peroxidase • 3 vials of dilution buffer • 3 vials of F.VII calibrator • 3 vials of washing solution • 3 vials of F.VII control	3 x 32 tests
00491	Asserachrom VIIa-AT Quantitative determination of factor VIIa Antithrombin complex by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of sample diluent • 3 vials of anti-AT-peroxidase • 3 vials of TMB • 3 vials of F.VIIa AT calibrator • 3 vials of F.VIIa AT control • 3 vials of AT peroxidase buffer • 1 vial of washing solution	3 x 32 tests
00956	Asserachrom X:Ag Quantitative determination of factor X:Ag by ELISA Method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-X:Ag-peroxidase • 3 vials of dilution buffer • 3 vials of F.X calibrator • 1 vial of washing solution • 3 vials of F.X control	3 x 32 tests
00280	Asserachrom VIII:Ag Quantitative determination of factor VIII:Ag by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-VIII:Ag-peroxidase • 3 vials of dilution buffer • 3 vials of F.VIII:Ag calibrator • 1 vial of washing solution • 3 vials of F.VIII:Ag control	3 x 32 tests
00264	Asserachrom sEPCR Quantitative determination of soluble Endothelial Protein C Receptor by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 1 vial of washing solution • 3 vials of dilution buffer • 3 vials of sEPCR calibrator • 3 vials of sEPCR control • 3 vials of anti-sEPCR-peroxidase	3 x 32 tests
00261	Asserachrom Total TFPI Quantitative determination of Total Tissue Factor Pathway Inhibitor (TFPI) by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of anti-Total TFPI-peroxidase • 1 vial of washing solution • 6 tablets of OPD • 3 vials of Total TFPI calibrator • 6 tablets of urea peroxide • 3 vials of Total TFPI control • 1 vial of dilution buffer	3 x 32 tests
00616*	Asserachrom TAFIa/TAFIai Quantitative determination of activated and/or inactivated TAFI by ELISA method <i>Research Product</i>	3 x 2 coated strips • 3 vials of TMB • 3 vials of anti-TAFIa/TAFIai Peroxidase • 3 vials of dilution buffer • 3 vials of TAFIa/TAFIai calibrator • 1 vial of washing solution • 3 vials of TAFIa/TAFIai control	3 x 32 tests
01004	Cy-Quant ELISA sCD146 Quantitative determination of soluble CD146 <i>Research Product</i>	3 x 2 coated strips • 6 tablets of OPD	3 x 32 tests
01514	Cy-Quant MV/TF Activity	1 vial of washing solution • 1 vial of non blocking TF Mab • 1 vial of anti-TF Mab • 2 vials of FVII + FX mix • 1 vial of dilution buffer • 1 vial of stop solution • 1 vial of chromogenic substrate • 2 vials of Calibrator • 2 vials of low activity control • 2 vials of high activity control	96 tests

*Availability depending upon the country.

Research kits

Cat. NR.	Product name	Product description	Pack.
Immuno turbidimetric methods			
00581	Liatest C4b-BP Quantitative determination of C4b-BP by immuno- turbidimetric method <i>Research Product</i>	6 vials of latex • 6 vials of buffer	6 x 1 mL
Flow Cytometry			
00452	Platelet PAIg Kit for Platelet Associated Immunoglobulin quantitation by flow cytometry <i>Research Product</i>	1 vial of diluent • 1 vial of calibrator • 1 vial of Mab • 1 vial of staining reagent • 1 vial of negative isotypic control • 1 vial of buffer	10 samples
00457	Platelet Calibrator Kit for customized platelet antigen quantitation by flow cytometry <i>Research Product</i>	1 vial of diluent • 1 vial of negative isotypic control IgG2b • 1 vial of negative isotypic control IgG1 • 1 vial of calibrator • 1 vial of negative isotypic control IgG2a • 1 vial of staining reagent	50 samples
00418	Platelet GP Screen Kit for customized platelet glycoprotein quantitation by flow cytometry <i>Research Product</i>	1 vial of diluent • 1 vial of Mab2 anti-GP Ia • 1 vial of calibrator • 1 vial of Mab2 anti-GP Ib • 1 vial of staining reagent • 1 vial of Mab2 anti-GP IIIa	10 samples
00420	Megamix Beads for cytometer settings in microparticles analysis <i>Research Product</i>	1 vial of beads	50 tests
01077	Megamix-Plus FSC Beads for FSC-optimized cytometer settings in microparticles analysis <i>Research Product</i>	1 vials of beads	50 tests
01078	Megamix-Plus SSC Beads for SSC-optimized cytometer settings in microparticles analysis <i>Research Product</i>	1 vial of beads	50 tests
01169	Mp-Count Beads Absolute counting microparticles by flow cytometry on PPP or purified microparticles <i>Research Product</i>	1 vial of 3 mL	100 tests

Purified proteins

Cat. NR.	Product name	Pack.
00461	Purified VWF	1 vial
00519	Purified Fibrinogen	1 vial
00557	Purified Prothrombin	1 vial
00896	Purified Thrombin (Human)	1 vial
00912	Purified Factor Xa	1 vial
00888	Purified AT III	1 vial
00828	Purified APC	1 vial
00392	Purified β 2 GlycoProtein I	1 vial
00964	Purified Bovine Thrombin	1 vial

Chromogenic substrates & Activators

Cat. NR.	Product name	Product description	Pack.
Chromogenic substrates			
00811	CBS 31.39 (Factor Xa)	Availability: 3 months maximum	1 vial
00873	CBS 34.47 (Thrombin)	Availability: 3 months maximum	1 vial
Activators			
00823	Ecarin	Availability: 3 months maximum	1 vial
00830	r-Hirudin	Availability: 3 months maximum	1 vial
00501	Ristocetin	Availability: 3 months maximum	1 x 5 mL

*Availability depending upon the country.

Monoclonal Antibodies for Flow Cytometry

Cat. NR.	Product name	Pack.
Anti-platelet markers		
01033	CD32, clone 2B2, purif.	0.1 mg
01031	CD36, clone 10.5, purif.	0.1 mg
01030	CD36, clone 10.5, FITC	100 tests
01032	CD36, clone 10.5, PE	100 tests
01025	CD41, clone PL2 49, purif.	0.1 mg
01024	CD41, clone PL2 49, FITC	100 tests
01026	CD41, clone PL2 49, PE	100 tests
01028	CD42b, clone ALMA 19, purif.	0.1 mg
01027	CD42b, clone ALMA 19, FITC	100 tests
01029	CD42b, clone ALMA 19, PE	100 tests
01017	CD61, clone 4F8, purif.	0.1 mg
01016	CD61, clone 4F8, FITC	100 tests
01022	CD62P, clone LYP20, purif.	0.1 mg
01021	CD62P, clone LYP20, FITC	100 tests
01023	CD62P, clone LYP20, PE	100 tests
01005	Fibrinogen, clone 9F9, FITC	100 tests
01083	GPVI, clone 1G5, purif.	0.1 mg
01084	GPVI, clone 1G5, PE	100 tests
Anti-endothelial cell markers		
01006	CD146, clone COM3D9, purif.	100 tests
01007	CD146, clone COM2F6, purif.	100 tests
01008	CD146, clone COM5G6, purif.	100 tests
01148	CD146, clone COM7A4, purif.	100 tests
01010	CD146, clone COM7A4, FITC	100 tests
01009	CD146, clone COM7A4, Biot.	100 tests
01149	CD146, clone S ENDO 1, purif.	100 tests
01013	CD146, clone S ENDO 1, FITC	100 tests
01015	CD146, clone S ENDO 1, PE	100 tests
01012	CD146, clone S ENDO 1, Biot.	100 tests
Isotypic Negative Controls		
01019	Ctl. neg. IgG1 purif. (2DNP2H11)	0.1 mg
01018	Ctl. neg. IgG1 FITC (2DNP2H11)	100 tests
01020	Ctl. neg. IgG1 PE (2DNP2H11)	100 tests
01038	Ctl. neg. IgG2a purif. (2DNP16C12)	0.1 mg
01037	Ctl. neg. IgG2a FITC (2DNP16C12)	100 tests
01039	Ctl. neg. IgG2a PE (2DNP16C12)	100 tests
01035	Ctl. neg. IgG2b purif. (2DNP14G5)	0.1 mg
01034	Ctl. neg. IgG2b FITC (2DNP14G5)	100 tests
01036	Ctl. neg. IgG2b PE (2DNP14G5)	100 tests

*Availability depending upon the country.

Thrombin Generation (CAT)

Cat. NR.	Product name	Product description	Pack.
86192	Thrombin Calibrator	20 vials	20 x 1 mL
86196	PRP Reagent	20 vials	20 x 1 mL
86222	MP Reagent	20 vials	20 x 1 mL
86193	PPP Reagent	20 vials	20 x 1 mL
86194	PPP Reagent Low	20 vials	20 x 1 mL
86195	PPP Reagent High	20 vials	20 x 1 mL
86197	FluCa kit	20 vials of Fluo buffer • 1 vial of Fluo substrate	20 x 1.6 mL 1 x 0.9 mL

Thrombin Generation (ST Genesis)

Cat. NR.	Product name	Product description	Pack.
01277	STG-Bleedscreen	3 vials of STG-Bleedscreen • 3 vials of STG-RefPlasma BLS • 3 vials of STG-QualiTest Norm BLS • 3 vials of STG-QualiTest Low BLS	3 x 4 x 1 mL
01279	STG-Thromboscreen	3 vials of STG-ThromboScreen -TM • 3 vials of STG-ThromboScreen+TM • 3 vials of STG-RefPlasma TS • 3 vials of STG-QualiTest High TS • 3 vials of STG-QualiTest Norm TS • 3 vials of STG-QualiTest Low TS	3 x 6 x 1 mL
01278	STG-Drugscreen	3 vials of STG-Drugscreen • 3 vials of STG-RefPlasma DS • 3 vials of STG-QualiTest Norm DS • 3 vials of STG-QualiTest Low DS	3 x 4 x 1 mL
01281	STG-Cal&Fluo	3 vials of STG-ThrombiCal • 3 vials of STG-FluoStart • 3 vials of STG-FluoSet	3 x 2 mL 3 x 1.5 mL 3 x 1.5 mL
01280	STG-Thrombiclean	6 vials of STG-ThrombiClean	6 x 2 mL

ExpertCor

Cat. NR.	Product name	Product description	Pack.
Correlation Plasma sets			
01494	ExpertCor Routine	30 vials for PT, APTT and fibrinogen • 5 vials for INR	35 x 1.5 mL
01307	ExpertCor DDI	20 vials	20 x 1 mL
01308	ExpertCor UFH	10 vials	10 x 1 mL
01309	ExpertCor LMWH	10 vials	10 x 1 mL

Heparin Induced Thrombocytopenia

Cat. NR.	Product name	Pack.
01482	5B9 Monoclonal Antibody	1 x 0.5 µg

*Availability depending upon the country.



INDEX

1		
1A ST Genesia fuse	17	
2		
2A ST Genesia fuse	17	
3		
3A ST Genesia fuse	17	
4		
4A ST Genesia fuse	17	
5		
5B9 Monoclonal Antibody	66	
5 x 20 T 10A fuses	6, 11	
5 x 20 T 1A fuses	6, 11	
5 x 20 T 2A fuses	6	
5 x 20 T 3.15A fuses	15	
5 x 20 T 315mA fuses	11	
5 x 20 T 4A fuses	11	
5 x 20 T 5A fuses	6	
5 x 20 T 6.3A fuses	6, 11, 13	
5 x 20 T 8A fuses	11	
6		
6.3 x 32 T 0.3A fuses	11	
6.3 x 32 T 10A fuses	6, 11	
6.3 x 32 T 15A fuses	6	
6.3 x 32 T 1A fuses	11	
6.3 x 32 T 2A fuses	6	
6.3 x 32 T 3A fuses	15	
6.3 x 32 T 4A fuses	11	
6.3 x 32 T 5A fuses	6	
6.3 x 32 T 6.25A fuses	6, 11	
6.3 x 32 T 8A fuses	11	
A		
ADP	50, 51	
Air filter	6, 11, 13	
Aluminium microcontainer reducer (microvolume)	7, 11	
Application notebook SStart Max	15	
Arachidonic acid	51	
Asserachrom D-Di	40	
Asserachrom Free Protein S	46	
Asserachrom HPIA	44	
Asserachrom HPIA-IgG	44	
Asserachrom IX:Ag	38	
Asserachrom PAI 1	42	
Asserachrom PF4	51	
Asserachrom Protein C	46	
Asserachrom sEPCR	47, 62	
Asserachrom TAFIa/TAFIai	42, 62	
Asserachrom Total Protein S	46	
Asserachrom Total TFPI	47, 62	
Asserachrom tPA	42	
Asserachrom VII:Ag	38, 62	
Asserachrom VIIa-AT	38, 62	
Asserachrom VIII:Ag	38, 62	
Asserachrom VWF:Ag	49	
Asserachrom VWF:CB	49	
Asserachrom VWF:FVIII B	38, 49	
Asserachrom Washing Solution	59	
Asserachrom X:Ag	38, 62	
Asserachrom β TG	51	
B		
Ball dispenser	15	
Ball extractor	7, 11, 13, 15	
Ball vial	15	
C		
Calibrated Automated Thrombogram	18	
C.K. Prest 2	32, 35	
C.K. Prest 5	32, 35	
CBS 31.39 (Factor Xa)	64	
CBS 34.47 (Thrombin)	64	
CD146, clone COM2F6, purif.	52, 65	
CD146, clone COM3D9, purif.	52, 65	
CD146, clone COM5G6, purif.	52, 65	
CD146, clone COM7A4, Biot.	52, 65	
CD146, clone COM7A4, FITC	52, 65	
CD146, clone COM7A4, purif.	52, 65	
CD146, clone S ENDO 1, Biot.	52, 65	
CD146, clone S ENDO 1, FITC	52, 65	
CD146, clone S ENDO 1, PE	52, 65	
CD146, clone S ENDO 1, purif.	52, 65	
CD32, clone 2B2, purif.	52, 65	
CD36, clone 10.5, FITC	52, 65	
CD36, clone 10.5, PE	52, 65	
CD36, clone 10.5, purif.	52, 65	
CD41, clone PL2 49, FITC	52, 65	
CD41, clone PL2 49, PE	52, 65	
CD41, clone PL2 49, purif.	52, 65	
CD42b, clone ALMA 19, FITC	52, 65	
CD42b, clone ALMA 19, PE	52, 65	
CD42b, clone ALMA 19, purif.	52, 65	
CD61, clone 4F8, FITC	52, 65	
CD61, clone 4F8, purif.	52, 65	
CD62P, clone LYP20, FITC	52, 65	
CD62P, clone LYP20, PE	52, 65	
CD62P, clone LYP20, purif.	52, 65	
Coag Control N+P	34, 35, 57	
Coag.One by Stago	22, 27	
Collagen	49, 51	
Colorimetric box filter	6, 11	
Combitips 2.5 mL (position 1 = 50 μ L)	15	
Complete dispensing tube assem.	18	
Connect.One by Stago	24	
Connect.One Gateway kit.	22	
Connected pipette	15	
Ctl. neg. IgG1 FITC (2DNP2H11)	53, 65	
Ctl. neg. IgG1 PE (2DNP2H11)	53, 65	
Ctl. neg. IgG1 purif. (2DNP2H11)	53, 65	
Ctl. neg. IgG2a FITC (2DNP16C12)	53, 65	
Ctl. neg. IgG2a PE (2DNP16C12)	53, 65	
Ctl. neg. IgG2a purif. (2DNP16C12)	53, 65	
Ctl. neg. IgG2b FITC (2DNP14G5)	53, 65	
Ctl. neg. IgG2b PE (2DNP14G5)	53, 65	
Ctl. neg. IgG2b purif. (2DNP14G5)	53, 65	
Cy-Quant ELISA sCD146	62	
Cy-Quant VASP/P2Y12	50	
Cy-Quant MV/TF Activity	50, 62	
D		
D-Di Test	40	
Dispensing tip	18	
E		
Ecarin	64	
Epinephrin	51	
Equipped needle	13, 17	
Etaloquick	34, 55	
ExpertCor	57, 66	
ExpertCor DDI	57, 66	
ExpertCor LMWH	57, 66	
ExpertCor Routine	57, 66	
ExpertCor UFH	57, 66	
F		
F.S. Test	40, 41, 59	
F.S. Test Unit	40	
Factor XIII Reagents	38	
FDP Plasma	40, 41	
FDP Plasma (buffer)	41	
FDP Plasma (latex)	41	
Fibrinogen, clone 9F9, FITC	52, 65	
Finntips 1.25 mL (position 1 = 25 μ L)	15	
Finntips 2.5 mL (position 1 = 50 μ L)	15	
FluCa kit	66	
G		
Glass tubes + stirring bars	19	
GPVI, clone 1G5, PE	52, 65	
GPVI, clone 1G5, purif.	52, 65	

H			
Halogen lamp	6, 11, 18		
Hamilton syringe and o-ring	6, 11		
I			
Immulon 2HB Plate	18		
K			
K-Assay Coagulation Control	36, 38		
K-Assay Factor XIII	36, 38		
K-Assay Factor XIII Calibrator	36, 38		
Kit of 3 filters ST Genesis	17		
L			
Liatest AT III	46		
Liatest C4b-BP	46, 63		
Liquid filter	6, 11, 13, 17		
M			
Megamix	50, 63		
Megamix-Plus FSC	50, 63		
Megamix-Plus SSC	50, 63		
Microcontainer standard rack	7		
Microcontainer	5, 10, 13		
Microcups reducer	7, 11, 13		
Mini USB Adaptor	15		
MP Reagent	66		
Mp-Count Beads	50, 63		
My Expert QC	26, 27		
My Qualiris QC	28, 29		
N			
Needle 1	6, 11		
Needle 2	6, 11		
Needle 3	6, 11		
Néoplastine CI Plus 2	32, 34		
Notebook STart Max support	15		
P			
PEG 25%	59		
Pipette thermal conductor grey	15		
Pipette thermal conductor blue	15		
Pipette thermal conductor black & grey	15		
Platelet Aggregometers TA-8V & TA-4V	19		
Platelet Calibrator	50, 63		
Platelet GP Screen	50, 63		
Platelet PAlg	50, 63		
PLT GP Receptors	50		
PLT VASP/P2Y12	50		
Pool Norm	32, 35, 47		
PPP Reagent	66		
PPP Reagent High	66		
PPP Reagent Low	66		
PRP Reagent	66		
PTT Automate 5	32, 35		
PTT-LA	32, 47		
Purified APC	64		
Purified AT III	64		
Purified Bovine Thrombin	64		
Purified Factor Xa	64		
Purified Fibrinogen	38, 64		
Purified Prothrombin	64		
Purified Thrombin (Human)	64		
Purified VWF	64		
Purified β 2 GlycoProtein I	64		
Q			
Qualiris Diagnostic Challenge	28, 58		
Qualiris QC D-Dimer	28, 58		
Qualiris QC Dabigatran	28, 58		
Qualiris QC DOAC Anti-Xa	28, 58		
Qualiris QC Experience	28, 58		
Qualiris QC FM-FDP	28, 58		
Qualiris QC Heparin HBPM/LMWH	28, 58		
Qualiris QC Heparin HNF/UFH	28, 58		
Qualiris QC Lupus Anticoagulant	28, 58		
Qualiris QC Premium Factors	28, 58		
Qualiris QC Premium S1	28, 58		
Qualiris QC Premium S2	28, 58		
Qualiris QC Premium Thrombophilia	28, 58		
QWERTY keyboard for STA Compact Max	11		
R			
r-Hirudin	64		
Rack for primary sample tubes ST Genesis	17		
Rack for secondary sample tubes ST Genesis	17		
Reagent tubing	18		
Reagents carousel n°1	13		
Reagents carousel n°2	13		
Red stirring magnet	5, 10		
Reducer for microvolume sample	13		
Reduction ring DIN 14 (18.5 mm)	13, 15		
Reduction ring DIN 18 (22.8 mm)	15		
Reducer for 5 mL vial	15		
Rising plate assay air filter	11		
Ristocetin	49, 51, 64		
S			
Sample carousel n°1	13		
Sample carousel n°2	13		
Sample rack	7		
Sample rack 11-13	7		
Sample rack holder	7		
Sample tray	7		
SCE upgrade to Coag.One server license	22		
SCE upgrade to Coag.One server license Workcell	22		
Screw microtube kit	7		
Soft stilet kit	6, 11		
SPA 20	32, 34		
ST Genesis system	16, 17		
STA Compact Max systems	9		
STA R Max systems	4, 5		
STA Satellite Max	12		
STA Satellite USB	13		
STA-Apixaban Calibrator	45, 55		
STA-Apixaban Control	45, 57		
STA-C.K. Prest 5	32, 35		
STA-CaCl ₂ 0.025 M	59		
STA-Cephascreen 10	32, 35		
STA-Cephascreen 4	32, 35		
STA-Cleaner Solution	5, 10, 13, 17, 59		
STA-Coag Control N+P	56		
STA-Control LA 1+2	56		
STA-Cuvettes	5, 10		
STA-D-Di Control	56		
STA-Dabigatran Calibrator	55		
STA-Dabigatran Control	57		
STA-Deficient II	36, 39		
STA-Deficient IX	37, 39		
STA-Deficient V	36, 39		
STA-Deficient VII	36, 39		
STA-Deficient VIII	37, 39		
STA-Deficient X	36, 39		
STA-Desorb U	5, 10, 13, 59		
STA-ECA II	44, 45		
STA-Edoxaban Calibrator	55		
STA-Edoxaban Control	57		
STA-FDP Calibrator	55		
STA-FDP Control	56		
STA-FM Calibrator	55		
STA-FM Control	56		
STA-Fondaparinux Calibrator	55		
STA-Fondaparinux Control	56		
STA-HBPM/LMWH Control	56		
STA-Hepanorm H	55		
STA-Hepanorm HBPM/LMWH	55		
STA-Heparin Control	56		
STA-ImmunoDef II	36, 39		
STA-ImmunoDef IX	37, 39		
STA-ImmunoDef VIII	37, 39		
STA-ImmunoDef XI	37, 39		
STA-ImmunoDef XII	37, 39		
STA-Liatest Control N+P	56		
STA-Liatest D-Di Plus	40, 41		
STA-Liatest FDP	41		
STA-Liatest FM	40, 41		
STA-Liatest Free Protein S 2	46, 48		
STA-Liatest Free Protein S 6	46, 48		
STA-Liatest VWF:Ag	49		
STA-Liquid Anti-Xa 4	43		
STA-Liquid Anti-Xa 8	43		
STA-Liquid Fib	33, 35		
STA-Maxi Reducer	5, 10, 13		
STA-Micro Reducer	5, 10		
STA-Microcups	5, 10, 13		
STA-Mini Reducer	5, 10, 13		
STA-Multi Hep Calibrator	55		
STA-NeoPTimal 10	32, 34		
STA-NeoPTimal 20	32, 34		
STA-NeoPTimal 5	32, 34		
STA-Néoplastine R 15	32, 34		
STA-Owren Koller	59		
STA-Procoag PPL	62		
STA-PTT Automate 5	32, 35		

STA-Quality HBPM/LMWH	56	STG-Cuvettes	17
STA-Quality HNF/UFH	56	STG-CuvetteWaste	17
STA-Reptilase	33, 35	STG-DrugScreen	54, 66
STA-Rivaroxaban Calibrator	55	STG-ThrombiClean	54, 66
STA-Rivaroxaban Control	57	STG-ThromboScreen	54, 66
STA-Routine QC 2 mL	56	STic Expert HIT 20	44
STA-Routine QC P Plus	56	STic Expert HIT 5	44
STA-Satellite Cuvettes	13	STic Expert HIT Control	57
STA-SPA+	32, 34	Suction head	6, 11
STA-Stachrom Antiplasmin	42	System Control N+P	57
STA-Stachrom AT III 3	46, 48	T	
STA-Stachrom AT III 6	46, 48	Teflon tips for syringe and o-ring	6, 11
STA-Stachrom Plasminogen	42	Test cards for D-Di Test kit	40, 59
STA-Stachrom Protein C	46, 48	Test cards for F.S. Test and FDP Plasma kits	40, 41, 59
STA-Stachrom TAFI	42, 62	Thermal paper	13, 15
STA-Staclot APC R	46	Thrombin Calibrator	66
STA-Staclot DRVV Confirm	47, 48	TRAP 6	51
STA-Staclot DRVV Screen 2	47, 48	Triangular stirring bar (6x12)	15
STA-Staclot DRVV Screen 5	47, 48	TriniCHROM Factor VIII:C	37, 38, 39
STA-Staclot Heparin 1	43	Tubing EV piercing	11
STA-Staclot Protein C 1	46, 48	U	
STA-Staclot Protein C 3	46, 48	Unicalibrator	55
STA-Staclot Protein S	46, 48	USB bar code handheld	15
STA-System Control N+P	56	V	
STA-Thrombin 10	33, 35	Visual Alarm Device	7
STA-Thrombin 2	33, 35	W	
STA-Unicalibrator	55	Wheel for tube rotation	6
STA-VWF:Ag Calibrator	55	White stirring magnet	5, 10, 13
STA-VWF:Rco	49		
Stachrom HCII	47, 62		
Stachrom Heparin	43		
Stachrom PAI	42		
Staclot LA	47, 48		
Staclot VIIa-rTF	38, 62		
STart Max	14		
Steel ball dispenser	15		
STG-BleedScreen	54, 66		
STG-Cal&Fluo	54, 66		

Asserachrom, Calibrated Automated Thrombogram (CAT), C.K. Prest, Cephascreen, Coag.One, Connect.One, CY-Quant, D-Di Test, Etaloquick, Hepanorm, K-Assay, My Expert QC, Néoplastine, Qualiris, ST Genesis, STA, STA Compact Max, STA-Liatest, STA-NeoPTimal, STA-R Evolution, STA R Max, STA Satellite, STA Workcell Max, STA Satellite Max, STG, Stachrom, Staclot, STart Max, STart 4, STic Expert, Thrombogram (CAT) are trademarks of the Stago Group.

The rights of the trademarks and logos used in this document belong to the Stago Group.
The use of these trademarks is not permitted without permission from the Stago Group.

NOTES

A series of horizontal dashed lines for writing notes.


Educational tools

Publishing

Practical Manual series

- One book per year since 2014
- Up to date clinical and lab practical information
- International scientific boards with renowned experts
- 6 issues currently available: Antiphospholipid Syndrome (**ref. 29289**) - Parenteral anticoagulant (**ref. 29618**) - Oral anticoagulant (**ref. 29691**) - Bleeding disorders (**ref. 300588**) - Platelet (**ref. 300726**) - Venous Thromboembolism (**ref. 301473**)

Available in English. Size: A5 - Pages: around 100



Quick guide

Hemostasis and Thrombosis essentials in your pocket: Screening assays, Anticoagulant therapy, DIC, Thrombophilia, D-Dimer assay, Lupus anticoagulants.

(**ref. 28171**)

Available in English, French, Italian, Portuguese and Spanish.



Practical guide to DOACs

Everything you need to know about direct oral anticoagulants (DOACs) in a posological sheet format: Indications - Dosing - Interactions - Contraindications - Pharmacokinetics - The drug's effect on hemostasis tests - Management of DOACs in various clinical situations - Stago solutions.

(**ref. 300345**)


Available in English, French, Italian, Portuguese and Spanish.



Reference ranges & Recommendations

Practical tools in a posological sheet format on various topics: Hemostasis and Pregnancy (**ref. 29211**) - Hemostasis and Pediatrics (**ref. 29202**) - Anticoagulation in children (**ref. 29201**) - ISTH diagnostic algorithm for overt DIC (**ref. 26905**).


Available in English, French, German, Italian, Portuguese, Spanish, depending on references.



History of Thrombosis/ History of Haemophilia

Twin posters depicting the history of 2 main coagulation defects: Thrombosis (**ref. 29628**) and Haemophilia (**ref. 300420**).

Available in Chinese, Dutch, English, French, German, Italian, Portuguese, Spanish, depending on posters.



All references mentioned are those of the English versions.



Diagnostica Stago S.A.S.
3 allée Thérèse - CS 10009
92600 Asnières sur Seine Cedex - France
+33 (0)1 46 88 20 20
+33 (0)1 47 91 08 91
webmaster@stago.com
www.stago.com